



Absorbable Beauty Non Surgical Face Lift Threading Mono 29G 50mm

Basic Information

- Place of Origin: china
- Brand Name: PDO PCL PLLA Threads
- Certification: ce
- Model Number: Mono 29G 50mm Pdo Thread
- Minimum Order Quantity: 10pcs
- Price: USD +0.5-6+pcs
- Packaging Details: 10pcs/1package
- Delivery Time: 3-5days
- Payment Terms: Western Union, T/T, D/P, D/A, L/C, MoneyGram
- Supply Ability: 10000pcs +pcs+week



Product Specification

- Size: 38mm 60mm 90mm, 19G 21G 23G 23G 26G 27G 29G 30G
- Color: Light Purple/white
- Feature: Easy To Operate, Two-year Expiry
- Packing: Sterile Packaging & Aluminium Foil Bag
- Properties: Injection & Puncture Instrument, Medical Polymer Materials & Products
- Highlight: **Absorbable non surgical face lift threading, non surgical face lift threading, beauty thread lift Mono 29G 50mm**



More Images



Our Product Introduction

for more products please visit us on hadermal-filler.com

Product Description

Medical Non-Surgical Absorbable Face Eyebrow Lifting Mono 29G 50mm Pdo Thread/ Mono Pdo Threads PDO PCL PLLA Threads

Which is better PDO or PCL thread?

PDO, or polydioxanone threads, dissolve after only six to nine months. Contrastingly, PLLA, or Poly L. Lactic Acid threads, dissolve after 12 to 18 months. PCL, or polycaprolactone threads, are the most durable of the three. What's the difference between PDO and PLLA threads?

The PLLA has been observed to be even more effective at simulating collagen production than PDO. This material has a lifetime of two years or more, and once dissolved, the PLLA breaks down into harmless ordinary substances like lactate, glucose, carbon dioxide and water. How long does a PDO thread lift last?

Is a thread lift worth it?

Efficacy. Thread lift procedures aren't as dramatically effective as facelifts, and studies on their long-term efficacy are lacking. Results from a thread lift last from 1 to 3 years. For best results, surgeons recommend combining a thread lift with other kinds of anti-aging procedures, such as ultrtherapy.

It almost sounds too good to be true, but an easy, minimally invasive treatment now makes it possible to lift and tighten skin, achieving results similar to a surgical facelift. A PDO Thread Lift leaves the face firmer, tighter, and more youthful in appearance, and it offers results that can last 18 months or more.


10pcs	1packge
shipping cost	25usd
diliverday	3-5days

PRODCTS DESCRIPTION





THREE TYPE OF NEEDLE



TYPES OF THREAD

PDO

PCL

NEW PRODUCT: DOUBLE THREAD

THE NEW TWO-WIRE DESIGN, THE LINE IS LONGER, AND THE PULLING EFFECT IS MORE OBVIOUS.

- **MONO**
- **Cog 1D**
- **Cog 2D**
- **Cog 3D**
- **Cog 4D**
- **Cog 6D**




||| SITE OF ACTION



COG WITH L CANNULA

COG3-1L(NOSE)	Needle		PDO
	Gauge	Length	Length
	19G	25mm	40mm
		50mm	70mm

COG4-1L	Needle		PDO
	Gauge	Length	Length
	19G	100mm	150mm
	21G	100mm	150mm

COG 3DL	Needle		PDO
	Gauge	Length	Length
	19G	100mm	150mm
	21G	100mm	150mm

COG 4DL	Needle		PDO
	Gauge	Length	Length
	18G	100mm	150mm
	19G	100mm	150mm
	21G	100mm	150mm

COG WITH W CANNULA

COG 4D W	Needle		PDO
	Gauge	Length	Length
	19G	100mm	150mm

COG 6D W	Needle		PDO
	Gauge	Length	Length
	18G	100mm	150mm

Tag Ring	Needle		PDO
	Gauge	Length	Length
	21G	38mm	50mm
	23G	38mm	50mm

Multi Thread	Needle		PDO
	Gauge	Length	Length
	21G	38mm	60mm
	23G	38mm	60mm

Mesh	Needle		PDO
	Gauge	Length	Length
	19G	38mm	50mm

Double Needle	Needle		PDO
	Gauge	Length	Length
	20G	100mm	210-560mm

Eye Thread	Needle		PDO
	Gauge	Length	Length
	30G	25mm	30mm
	30G	38mm	60mm

Nose Thread	Needle		PDO
	Gauge	Length	Length
	19G	38mm	50mm
	21G	38mm	50mm

Fish Bone Thread	Needle		PDO
	Gauge	Length	Length
	18G	100mm	160mm
	19G	100mm	160mm



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA