



18G 19G Micro Cannula Needle 14-30g For Derma Filler Injection

Our Product Introduction

Basic Information

- Place of Origin: korea
- Brand Name: Hypodermic
- Certification: ce
- Model Number: 14g 16g 18G 19G 21G 23G 25g 27g 29 30G 31g 32g 34g
- Minimum Order Quantity: 1 box
- Price: usd 40
- Packaging Details: 50pcs/box
- Delivery Time: 3-5day
- Payment Terms: Western Union, MoneyGram, T/T,credit card ,paypal
- Supply Ability: 100000box



Product Specification

- Product Name: What Is A Micro-Cannula? A Micro-cannula Is A Very Thin, Tiny, Flexible Hollow Tube With A Smooth Round Tip. The Sharp End Is Gone! Instead Of Many Needle Sticks, We Only Need A Very Few To Provide Entry Points For The Cannula. That Means Fewer Needle Mark
- Material: Stainless Steel
- Tip Type: Blunt Type
- Specification: 14-30g
- Packing: 50 Pcs/box
- Highlight: **18G 19G Micro Cannula Needle, Micro Cannula Needle 30g, Blunt derma filler injection**



More Images



for more products please visit us on hadermal-filler.com

Product Description

18G 19G 21G 23G 30G Hypodermic disposable blunt micro cannula needle derma filler injection


What Is a Micro-Cannula?

A micro-cannula is a very thin, tiny, flexible hollow tube with a smooth round tip. The sharp end is gone! Instead of many needle sticks, we only need a very few to provide entry points for the cannula. That means fewer needle marks, less pain and far less probability of bruising.

What Are the Benefits of a Micro-Cannula over a Needle?

- Reduced pain
- Less bleeding
- Less bruising
- Less tissue trauma or swelling

PRODUCT SIZE



Micro-cannula injection requires a single entry point for a large area instead of many as occurs with a traditional needle. Because the tip of the cannula is smooth and rounded, small vessels are simply pushed out of the way instead of being cut – so instead of pain you'll most frequently experience little more than pressure.

spec

Product Name	disposable micro cannula needle
material	stainless steel
tip type	blunt
spec	14g-34g
application	medical beauty clinic
size	38mm 50mm 70mm
Quality Certification:	ce
Product Name:	Disposable Plastic Micro Cannula Blunt Tip Needles
Properties	Injection & Puncture Instrument
Application:	Anti-Puffiness, Skin Tightening

How Does the Micro-Cannula Get into the Skin?

An entry point is made with a traditional needle that is slightly larger than the microcannula. If desired, Dr. Fitzgerald will use a small amount of lidocaine to numb that area. Then the microcannula slides under the skin to the areas to be treated with filler.

Mode	Specification		Color	Shape
	Metric	Imperial		
0101	0.4*30	27G	Grey	Round Tip
0102	0.4*40	27G	Grey	Round Tip
0103	0.5*38	25G	Orange	Round Tip
0104	0.5*50	25G	Orange	Round Tip
0105	0.6*50	23G	Dark Blue	Round Tip
0106	0.6*70	23G	Dark Blue	Round Tip
0107	0.7*50	22G	Black	Round Tip
0108	0.8*50	21G	Dark Green	Round Tip

Note: customization is acceptable

How Many Entry Points Will I Require?

For an upper and lower lip, only two – one at each end of the lip area. These same entry points can be utilized for filler in additional areas around the mouth and even into the lower cheek area. For a full lower face rejuvenation (cheeks, lip, mouth area, chin) you will likely have 4 points of entry.

Needles push the filler in front of the needle tip through an open end. Micro-cannulas have a small opening on the side and a closed end. This opening provides a smoother delivery of the filler so no matter what the pressure from the injector, there is a smoother flow of material.

PRODUCT DOSAGE





shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA