

Mesotherapy RF Beauty Machine 1.5" OLED Meso Injector Gun 5 Level

Our Product Introduction

Basic Information

- Place of Origin: china
- Brand Name: meso gun
- Certification: ce
- Model Number: meso
- Minimum Order Quantity: one box
- Price: USD+80-95+BOX
- Packaging Details: ONE BOX ONE PCS
- Delivery Time: 3-5WORK DAY
- Payment Terms: L/C, D/A, D/P, T/T, Western Union
- Supply Ability: 10000pcs +pcs+week



Product Specification

- Application: For Commercial
- Product Name: Mesotherapy Gun
- Feature: Wrinkle Remover, Skin Tightening, Anti-aging, Skin Revitalizer, Pigment Removal, Skin Rejuvenation, Wrinkle Remover
- Injection Setting: 5 Level (Continuous), 6 Level (Dose)
- Power: Battery Embedded, Lithium Polymer Battery
- Application: For Commercial
- Selling Units: Single Item
- Highlight: **Mesotherapy RF Beauty Machine, meso injector gun 5 level, 1.5" OLED meso injector gun**



More Images



for more products please visit us on hadermal-filler.com

Product Description

mesotherapy inmeso gun vacuum injection gun water mesotherapy tion for skin lift/mesotherapy gun RF Beauty MACHINE

Description:

Injection System

Power: Battery embedded, Lithium Polymer Battery (7.4V, 850mAh)

Charger: 110 ~ 240V 50/60Hz (Adapter Input), DC 8.4V 500mA Max (Output)

Duration time: The use of negative pressure up to 2,000 times.

Power consumption: 4W

Material: ABS (Main machine, Negative pressure handle, Charging base), PC (Syringe holder, Filter),

Nickel-plated (Syringe piston), Silicone tube (Negative Pressure Tube)

Display: 1.5 inch color OLED

Power Operation: Switch button or No automatic shutdown in 5 minutes.

Mode: Continuous, Dose

Injection Setting: 5 level (continuous), 6 level (dose)

Weihtg: Main body 220g (Include syringe holder), Negative pressure handle 100g

Package Included:

Negative Pressure Handle x 1

Main body x 1

Syringe Piston : 5cc

Charging Base x 1

Charger x 1

Syringe Holder : 1 ~ 3cc, 5 ~ 10cc

Filter Tip x 1

Negative Pressure Tube x 2

Negative Pressure Needle 5 pins x 1

Charge Converter x 1

Use manual x 1

Water Syringe Injector ---- The injector gun was developed for water mesotherapy - hyaluronic acid is injected into the skin to solve any problems caused by the lack of moisture.

Function---- It makes the skin smooth, soft, radiant and hydrated after the injection. In addition, it is also suitable to inject various skin care products

Color OLED display---- It can provide a complete vision from all angles. Removable docking vacuum handle makes the water easier injection more convenient


Effective Beauty Device---- The recommended injection technique is the micro-puncture technique, in which very small amounts of hyaluronic acid, VC or other skin boosters are injected into the deep dermis, about 1-2 cm apart

2 Types of Storage Keys---- M1 for water mode, which can shrink pores and block moisture, M2 for hyaluronic acid mode, which can help the skin to remove cosmetic residues



Two head mix fillers

Could mix two style fillers with this,
Convenient and easy to operate



Hydrating

Add moisture to the skin,
let the skin become moist

Water and oil balance

The skin does not dry
the good health of the state

Skin Tightening

Let loose skin rejuvenation,
compact

Bright color of skin

Smooth and fair skin






Before



After



Before



After



Before



After





shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA