



Corsslink Hyaluronic Acid Dermal Filler Korea Gel 24mg/ml

Our Product Introduction

for more products please visit us on hadermal-filler.com

Basic Information

- Place of Origin: china
- Brand Name: comefad
- Certification: ce
- Model Number: 1ml 2ml 3ml 5ml 10ml 20ml
- Minimum Order Quantity: 1ml
- Price: usd 14
- Packaging Details: one box 1 pcs
- Delivery Time: 3-5day
- Payment Terms: T/T, Western Union, MoneyGram,credit card ,paypal
- Supply Ability: 100000pcs/day



Product Specification

- Product Name: Hyaluronic Acid Gel,Hyaluronic Acid Dermal Filler,hyaluronic Acid Dermal Injection Filler 10 MI For Skin Tightening,Lip Dermal Filler HA Injection Cosfiller Brand Hyaluronic Acid Gel
- Function: Remove Wrinkle,anti-aging,breast Enhancement
- Material: Cross-linked HA,Hyaluronic Acid Filler,Non Animal Cross-linked Hyaluronic Acid
- Application: Breast And Buttock Enlargement,Body Filler,Lip Enhancement,Anti-Wrinkle
- Duration: 6-18 Months
- Shelf Life: 2 Years
- HA Concentration: 24mg/ml
- Packing: Box
- Highlight: **Corsslink Hyaluronic Acid Dermal Filler, Hyaluronic Acid Dermal Filler 24mg/ml**



More Images



Product Description

Corsslink Hyaluronic Acid Dermal Filler Korea Gel

1ml 2ml 3ml 5ml 10ml 20ml Cross Linked Injectable Hyaluronic Acid Facial Dermal Filler Lip Buttock

Dermal fillers are gel-like substances that are injected beneath the skin to restore lost volume, smooth lines and soften creases, or enhance facial contours. More than 1 million men and women annually choose this popular facial rejuvenation treatment, which can be a cost-effective way to look younger without surgery or downtime.



spec

Product name	Hyaluronic Acid Gel, Hyaluronic Acid Dermal Filler, hyaluronic acid dermal injection filler 10 ml for skin tightening, Filler for skin, Lip dermal filler HA injection Cosfiller Brand Hyaluronic acid Gel
Function	Cross Linked, safety, easy to use
Cross-linked agent	BDDE
Packing	Blister Package, 1 syringe/box, 1pc/box, One Prefilled Syringe with Two Needles.
Material	Pure Hyaluronic Acid, Cross-linked Sodium Hyaluronate, Cross-linked HA, Hyaluronic acid filler, Non Animal Cross-linked Hyaluronic Acid
Application	breast and buttock enlargement, Body Filler, Lip enhancement, Neck, Anti-Wrinkle
Type	fine,derm.deep.ultra deep
Ha concentration	24mg/ml, 20mg--30mg/ml
Particle size	0.28~0.5mm, 0.1mm-2.0mm
Duration	6~24 months, 6-18 months, 6-18months, 8~12 Months

What is injectable hyaluronic acid?

Injectable hyaluronic acid (HA) is a type of temporary dermal filler. Hyaluronic acid is found naturally through out the body, with the highest concentrations in the joints, eyes, and skin. In the skin hyaluronic acid attracts water, which allows the skin to become more hydrated. As we age, the fat, muscles, bone, and skin in our face begins to thin. This loss of volume leads to either a sunken or sagging appearance of the face, fine lines, wrinkles, folds, and thin lips. Injectable HA is used to reduce the appearance of fine lines and wrinkles, facial folds, and to create structure, framework, and volume to the face and lips. The effects of injectable hyaluronic acid are seen immediately.

Hyaluronic Acid Dermal Injection Filler				
Type	Fine	Derm	Deep	Subskin
Concentration of HA	24mg/ml	24mg/ml	24mg/ml	24mg/ml

Approximate Number of gel partide /ml	200,000	100,000	10,000	5000
Volume of Syringe	1ml/2ml	1ml/2ml	1ml/2ml	10ml/20ml
Partide Size(mm)	0.1-0.15mm	0.15-0.28mm	0.28-0.5mm	0.5-1.25mm
Recommended treatment	Thin superficial lines, such as worry lines, periobital lines, perioral lines	Moderate wrinkles, such as glabellar, oral commissures. lips: fuleness, pouting and vermilion border	Deep facial wrinkles and folds, like the nasolabial folds, shaping facial contours, eg cheeks fullness, chin and lips augmentation	The local depression, soft tissure filler, also can be used in rhinoplasty For breast and buttock enhancement
Application Area	Facial skin	Facial skin	Facial skin	Facial skin
Lasting Duration	6~12 months	6~12 months	8~12 months	10~18 months
Injection Depth	Upper part of dermis	Middle part of dermis	Deep layer of dermis and / or surface layer of subcuits	Deep dermis or subcutaneous tissue shallow
Needle Specification	BD needs 30g	BD needle 27g	BD needle 26g	BD needle 23g
Ingredient	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid
Storage	Room Temperature			
Origin	Non Animal cross linke Hyaluronic Acid			

Fine Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.1-0.15mm

Volume: 1ml/2ML

Needle: 30G

Duration: 9-15 Months

Indication: Thin superficial lines, such as shallow lines removal, Crow's feet, perioral lines, smile lines,etc

Derm Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.15-0.28mm

Volume: 1ml/2ML

Needle: 27G

Duration: 9-15 Months

Indication: Moderate wrinkles, such as glabellar, oral commissures, nasolabial folds

Lips: fullness, pouting and vermilion border

Deep Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.28-0.5mm

Volume: 1ml/2ML

Needle: 26G

Duration: 9-15 Months

Indication: Deep facial wrinkles and folds, restore cheeks volume, reshaping chin and jawline correction

Sub-skin Line Dermal Filler 10ml/20ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.5-1.25mm

Volume: 10ml/20ml

Needle: 23G

Duration: 12-18 Months

Indication: The local depression, soft tissue filler, also can be used in rhinoplasty

For breast and buttock enhancement

Specification

Used to treat signs of aging

	Fine Lines	Derm	Derm Deep	Ultra-Deep
Partical Size	0.10-0.15mm	0.15-0.28mm	0.28-1.25mm	1.25-2.0mm
Partical No pergel	300,000	100,000	10,000	500
Concentration	24mg/ml	24mg/ml	24mg/ml	24mg/ml
Volume	1ml&2ml	1ml&2ml	1ml&2ml	10ml
Neekle Size	30G	30G	27G	/
Duration	6-9 months	6-12 months	6-12 months	9-18 months
Indications	Correction of thin superficiaines such as worry lines,periorbital lines,periorlnes	Correction of moderate wrinkles, such as glabellar, lips fullness,pouting vermilion border	Correction of deep facial wrinkles and fold,shaping facial contours,cheeks fullness,chin.	Rhinoplasty,breast and buttock augmentation

How does the treatment work?

The treatment can be done in 15 minutes but it can also take longer, depending on the area to be treated.

The professional in charge of your treatment injects the product using an ultra-thin needle under the surface of the skin. For fine lines, the product is injected at the surface of the skin. For marked wrinkles, the injection is made at a deeper level.

Please note that the injection of Juvéderm Ultra requires a medical gesture that can only be carried out by a doctor (doctor specialised in dermoesthetics, dermatologist or plastic surgeon).



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA