

PDO PCL PLLA Threads for Face Lifting 4D 6D 8D 18g 19g 21g 100mm L W Blunt Type

Our Product Introduction

Basic Information

- Place of Origin: korea
- Brand Name: comefad
- Certification: ce
- Model Number: pdo pcl plla
- Minimum Order Quantity: 10 pcs
- Price: 0.5-5 dollar
- Packaging Details: Sterile Aluminum packaging
- Delivery Time: 3-5day
- Payment Terms: T/T, Western Union, MoneyGram, paypal, credit card
- Supply Ability: 1000000 bag/day



Product Specification

- Size: 18g 19g 21g 22g 23g 25g 26g 27g 28g 30g
- Shelf Life: 2 Years
- Thread Length: 25mm To 300mm
- Needle Type: Sharp/Blunt Cannula
- Packing: Sterile Aluminum Packaging
- Material: PDO/PCL/PLLA
- Used For: Aesthetics Beauty Clinic Spa
- Function: Face Lift, Body Lift, Skin Tightening
- Packing Size: 20pcs/bag
- Sterilization: Ethylene Oxide
- Usage: Face/Neck/Body
- Keywords: Pdo Threads
- Thread Type: Barbed/Cogged/Screw
- Thread Size: 29G/30G/31G
- Stock: Yes




More Images



for more products please visit us on hadermal-filler.com

Product Description

Lifting Threads PDO PCL PLLA 4D 6D 8D 18g 19g 21g 100mm L W Blunt Type for face lifting



What are PDO Threads? PDO, which stands for polydioxanone, is an absorbable polymer that is flexible and durable.. The threads can be thought of as sutures, and there are three different types: Mono threads, smooth, and cog. Mono threads are smooth and singular sutures that provide a small lift. Cog threads have more of a barbed texture, provide a bit more of a lift. Screw threads have multiple threads intertwined and provide more volume to the treated area. The best part about these threads – they will stimulate natural collagen production. Here at Regenesis, we offer Mint PDO threads, as well as Nova Threads.

SMOOTH

Image	Needle		PDO		Image	Needle		PDO		
	Gauge	Length	Length			Gauge	Length	Length		
	26G	60mm	90mm		25G	90mm	135mm	26G	38mm	50mm
		90mm	135mm			60mm	90mm			
	27G	38mm	50mm		27G	38mm	50mm	29G	38mm	50mm
		50mm	70mm			50mm	70mm			
	29G	38mm	50mm		29G	38mm	50mm	30G	25mm	30mm
		50mm	70mm			60mm	90mm			
	60mm	90mm		50mm	70mm		38mm	50mm		
MONO										
	30G	25mm	30mm		26G	90mm	135mm	27G	38mm	50mm
		38mm	50mm			50mm	70mm			
						60mm	90mm			
MONO BLUNT(EYE)										

COG

Image	Needle		PDO		Image	Needle		PDO		
	Gauge	Length	Length			Gauge	Length	Length		
	21G	60mm	110mm		19G	100mm	150mm	21G	100mm	150mm
		100mm	150mm			100mm	150mm			
	19G	100mm	150mm		18G	100mm	150mm	19G	100mm	150mm
		60mm	110mm			21G	100mm		150mm	
		100mm	150mm							
COG 3D										


COG WITH L CANNULA

Image	Needle		PDO		Image	Needle		PDO		
	Gauge	Length	Length			Gauge	Length	Length		
	19G	25mm	40mm		18G	100mm	150mm	19G	100mm	150mm
		50mm	70mm			100mm	150mm			
	19G	100mm	150mm		19G	100mm	150mm	21G	100mm	150mm
		60mm	110mm			21G	100mm		150mm	
		100mm	150mm							
COG 3D L										

COG WITH W CANNULA


Needle	PDO	
	Gauge	Length
COG 4D W 	19G	100mm
COG 6D W 	18G	100mm


Needle	PDO	
	Gauge	Length
Tag Ring 	21G	38mm
	23G	38mm


Needle	PDO	
	Gauge	Length
Multi Thread 	21G	38mm
	23G	38mm

Needle	PDO	
	Gauge	Length
Mesh 	19G	38mm

Needle	PDO	
	Gauge	Length
Double needle 	20G	100mm
	21G	160mm


Needle	PDO	
	Gauge	Length
Eye thread 	30G	25mm
	30G	38mm

Needle	PDO	
	Gauge	Length
Nose thread 	19G	38mm
	21G	38mm

Needle	PDO	
	Gauge	Length
Fish bone thread 	18G	100mm
	19G	100mm

how long do the pdo thread last?

The optimal candidate is a person in the 30 to 60 year old age group who wants to tighten or lift certain areas of skin but who isn't ready for a facelift surgery. Most people like to target the nasolabial folds that go on either side of the nose to the sides of the mouth, and the marionette lines that go from the corners of the mouth down to the side of the chin. But many other areas can be targeted including the cheeks, brow, jowls, neck and underneath the chin. People older than 60 can be treated but typically don't see results as good as those younger. Frequently, more senior adults and those seeking more dramatic results would benefit more from a surgical facelift



spec	
Size	18g 19g 21g 22g 23g 25g 27g 28g 29g 30g
Shelf Life:	2 years
Quality Certification:	ce
Application:	skin lift, anti age, wrinkle remove, skin tightening

Thread Length:	25mm -300mm
needle type	sharp /blunt needle
Material	pdo pcl pla
product name	pdo thread lift
packing	Aluminum packaging
Properties:	Medical Adhesive & Suture Material

PACKAGE



✓ MOQ : 1 Pack

✓ 2 Years Expiry Period



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA