



## 14g 16g 17g 18g 19g 38mm 50mm Micro Blunt Cannula Syringe Needle 23G 25g

### Our Product Introduction

for more products please visit us on [hadermal-filler.com](http://hadermal-filler.com)

#### Basic Information

- Place of Origin: China
- Brand Name: Micro Cannula Needle
- Certification: CE
- Model Number: 18G 21G 23G 25G 27G 30G
- Minimum Order Quantity: 10
- Price: 0.7
- Packaging Details: Sterile package
- Delivery Time: 1-3Days
- Payment Terms: T/T, Western Union
- Supply Ability: 80000



#### Product Specification

- Product Name: Micro Cannula Needle
- Usage: Single-Use
- Type: Blunt Tip
- Size: 25mm 38mm 50mm 70mm
- Disinfecting Type: Far Infrared
- Material: Stainless Steel
- Function: Injection & Puncture Instrument
- Warranty: 2 Years
- Spe: 14g 16g 17g 18g 19g



#### More Images



## Product Description

### Disposable Stainless Blunt Tip Needle Micro Cannula Needle 18G 21G 27G

#### What Is A Blunt Cannula?

Cannulae are blunt, which mean that they are less likely to pierce a through a vein or artery. This results in less swelling, bruising and intravascular injection of filler. The practitioner has more feel about what structures are being passed by the tip of the cannula.

#### What Is A Blunt Plastic Cannula Used For?

The BD blunt plastic cannula is a truly universal cannula designed to be used with all types of split septum IV injection sites and vials for needle-free access.


#### How Do You Use A Blunt Tip Cannula?

A needle is required to make an entry point for the cannula that just passes through the skin. The cannula is then used to pass through the entry point to the areas beneath the skin. The blunt-tipped cannula can re-enter the same entry point again without having to hurt the patient with another needle puncture.

Name	blunt cannula needle
size	15g 16g 18g 19g
size2	38mm 50mm 60mm 90mm



## Product specifications



The diagram shows several needles of different gauges and lengths. Callouts indicate specific features: 'Micro tip' points to the needle tip, 'Scale mark' points to a mark on the needle shaft, and a comparison of 'Regular wall' and 'Ultra thin wall' is shown in a circular inset.

Site	Length(mm)	Site	Length(mm)
18G	50	25G	38
	70		50
	100		70
21G	50	27G	38
	70		50
22G	50		70
	70	30G	38
23G	38		40
	50		50
	70		70

## Product features

**A/**

Flexible cannulas for superficial & stiffer cannulas for deep injections.

Blunt tips slide under the skin without causing any trauma.


**/B**




Specifically designed for atraumatic intradermal injections of fluids.




## Product advantages



Blunt-tip





Flexible



No pin

## Product shows



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA