₫

more products please visit us on hadermal-filler.com

# shijiazhuang tiantian technology co.,ltd

### Wholesale Stainless Steel Dermal Filler Facial Filling 22g 23G 25g 38mm 50mm Flexible Fine Blunt Micro Tip Cannula Needl

### Basic Information

• Place of Origin:

• Brand Name: Micro Cannula Needle

Certification: CF

18G 21G 23G 25G 27G 30G Model Number:

 Minimum Order Quantity: 10 • Price: 0.7

 Packaging Details: Sterile package Delivery Time: 1-3Days • Payment Terms: T/T, Western Union

Supply Ability: 80000



### **Product Specification**

• Product Name: Micro Cannula Needle

• Usage: Single-Use • Type: Blunt Tip

Size: 25mm 38mm 50mm 70mm

Disinfecting Type: Far Infrared Material: Stainless Steel

• Function: Injection & Puncture Instrument

Warranty: 2 Years

14g 16g 17g 18g 19g Spe:



### More Images











### Product Description

Disposable Stainless Blunt Tip Needle Micro Cannula Needle 18G 21G 27G

What Is A Blunt Cannula?

Cannulae are blunt, which mean that they are less likely to pierce a through a vein or artery. This results in less swelling, bruising and intravascular injection of filler. The practitioner has more feel about what structures are being passed by the tip of the cannula.

What Is A Blunt Plastic Cannula Used For?

The BD blunt plastic cannula is a truly universal cannula designed to be used with all types of split septum IV injection sites and vials for needle-


How Do You Use A Blunt Tip Cannula?

A needle is required to make an entry point for the cannula that just passes through the skin. The cannula is then used to pass through the entry point to the areas beneath the skin. The blunt-tipped cannula can re-enter the same entry point again without having to hurt the patient with another needle puncture.


Name	blunt cannula needle
size	15g 16g 18g 19g
size2	38mm 50mm 60mm 90mm



# Product specifications





# Product features



# Product advantages 1 Virtually no bruising 2 Lessen pain considerably New collagen Single point of entry for wider coverage Minimal residence time Rewible

# Product shows





shijiazhuang tiantian technology co.,ltd



● 008619931162089 info@tiantianbike.com hadermal-filler.com



NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA