



shijiazhuang tiantian technology co.,ltd  
hadermal-filler.com



## Buttock Augmentation 20ml 50ml 100ml Derm Plus Hyaluronic Acid Filler For Breast Enlargement Breast Injection Price

### Our Product Introduction

for more products please visit us on [hadermal-filler.com](http://hadermal-filler.com)

#### Basic Information

- Place of Origin: China
- Brand Name: WYKJ
- Certification: CE
- Model Number: Dermal Filler
- Minimum Order Quantity: 1 box
- Price: 20
- Packaging Details: box
- Delivery Time: 5-8days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 50000box/month



#### Product Specification

- Product Name: Hyaluronic Acid Filler
- Type: Fine/Derm/Deep/Ultra Deep
- Grade: 1ml 2ml 3ml 5ml 10ml 20ml
- Application: Facial Contour And Body
- Duration: 6-18 Months
- Shelf Life: 2 Years
- HA Concentration: 24mg/ml
- Packing: Box



#### More Images



## Product Description

10ml Hyaluronic Acid Enlargement Products Breast Buttock Dermal Filler Injection

### 2ml Anti Wrinkle Hyaluronic Acid Facial Lip Dermal Filler

#### What is injectable hyaluronic acid?

Injectable hyaluronic acid (HA) is a type of temporary dermal filler. Hyaluronic acid is found naturally through out the body, with the highest concentrations in the joints, eyes, and skin. In the skin hyaluronic acid attracts water, which allows the skin to become more hydrated. As we age, the fat, muscles, bone, and skin in our face begins to thin. This loss of volume leads to either a sunken or sagging appearance of the face, fine lines, wrinkles, folds, and thin lips. Injectable HA is used to reduce the appearance of fine lines and wrinkles, facial folds, and to create structure, framework, and volume to the face and lips. The effects of injectable hyaluronic acid are seen immediately.

Hyaluronic Acid Dermal Injection Filler				
Type	Fine	Derm	Deep	Subskin
Concentration of HA	24mg/ml	24mg/ml	24mg/ml	24mg/ml
Approximate Number of gel partide /ml	200,000	100,000	10,000	5000
Volume of Syringe	1ml/2ml	1ml/2ml	1ml/2ml	10ml/20ml
Partide Size(mm)	0.1-0.15mm	0.15-0.28mm	0.28-0.5mm	0.5-1.25mm
Recommend ed treatment	Thin superficial lines, such as worry lines, periorbital lines, perioral lines	Moderate wrinkles, such as glabellar, oral commissures. lips: fullness, pouting and vermilion border	Deep facial wrinkles and folds, like the nasolabial folds, shaping facial contours, eg cheeks fullness, chin and lips augmentation	The local depression, soft tissue filler, also can be used in rhinoplasty For breast and buttock enhancement
Application Area	Facial skin	Facial skin	Facial skin	Facial skin
Lasting Duration	6~12 months	6~12 months	8~12 months	10~18 months
Injection Depth	Upper part of dermis	Middle part of dermis	Deep layer of dermis and / or surface layer of subcuits	Deep dermis or subcutaneous tissue shallow
Needle Specification	BD needles 30g	BD needle 27g	BD needle 26g	BD needle 23g
Ingredient	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid	Non Animal Cross Link Hyaluronic Acid
Storage	Room Temperature			
Origin	Non Animal cross linke Hyaluronic Acid			

#### Fine Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.1-0.15mm

Volume: 1ml/2ML

Needle: 30G

Duration: 9-15 Months

Indication: Thin superficial lines, such as shallow lines removal, Crow's feet, perioral lines, smile lines,etc

#### Derm Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.15-0.28mm

Volume: 1ml/2ML

Needle: 27G

Duration: 9-15 Months

Indication: Moderate wrinkles, such as glabellar, oral commissures, nasolabial folds

Lips: fullness, pouting and vermilion border

#### Deep Line Dermal Filler 1ml/2ml

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.28-0.5mm

Volume: 1ml/2ML

Needle: 26G

Duration: 9-15 Months

Indication: Deep facial wrinkles and folds, restore cheeks volume, reshaping chin and jawline correction

### **Sub-skin Line Dermal Filler 10ml/20ml**

Ingredient: Crosslinked hyaluronic acid gel

HA Concentration: 24mg/ml

Particle size: 0.5-1.25mm

Volume: 10ml/20ml

Needle: 23G

Duration: 12-18 Months

Indication: The local depression, soft tissue filler, also can be used in rhinoplasty

For breast and buttock enhancement








FILLING  
FINE LINES

FIRMNESS  
ELASTICITY

BOUNCING  
TIGHTENING

## PARAMETER

				
	<b>Fine Line</b>	<b>Derm Line</b>	<b>Deep Line</b>	<b>Ultra Line</b>
<b>Recommended indications</b>	For correction of thin lines across fore-head and around eyes and mouse	For treatment of wrinkles and lip enhancement	For creating fuller lips and shaping facial contours, such as cheeks and chin	For breast enlargement and buttock augmentation
<b>Concentration</b>	20mg/ml stabilized hyaluronic acid	20mg/ml stabilized hyaluronic acid	20mg/ml stabilized hyaluronic acid	20mg/ml stabilized hyaluronic acid
<b>Particle size</b>	0.10-0.15(ml)	0.15-0.28(ml)	0.28-0.5(ml)	0.5-1.25(ml)
<b>Approximate number of gel particles</b>	150.000ml	100.000ml	100.000ml	5.000ml
<b>Needles</b>	30g*1/2*	27g*1/2*	27g*1/2*	23g*1/2*
<b>Duration</b>	3-9months	6-12months	6-12months	6-12months
<b>Use for area</b>				

## Injectable Filler Treatment Areas

Cheek Volume Loss

Acne Scars

Smile Lines

Marionette Lines

Lip Volumizing

Lip Lines

Chin & Jawline Contouring



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA

