



## 2 In 1 Hyaluron Pen 0.3ml/ 0.5ml Hyaluronic Acid Pen

Our Product Introduction

for more products please visit us on [hadermal-filler.com](http://hadermal-filler.com)

### Basic Information

- Place of Origin: china
- Brand Name: hyaluron pen
- Certification: ce
- Model Number: 0.3ml/0.5ml
- Minimum Order Quantity: 1
- Price: usd 30
- Packaging Details: 1 0.3ml ampoule,1 0.5ml ampoule,1 pen,1 box
- Delivery Time: 3-5day
- Payment Terms: Western Union, MoneyGram, T/T,paypal,credit card
- Supply Ability: 100000 pcs



### Product Specification

- Packing Size: 22.0 \* 12.0 \* 7.0 ( Cm )
- Applications: Wrinkle Remover And Lip Lifting
- Content: 1 0.3ml Ampoule,1 0.5ml Ampoule,1 pen,1 Box
- Technology: Needle Free High Pressure Atomization
- Highlight: **0.5ml Hyaluronic Acid Pen,**  
**2 in 1 Hyaluronic Acid Pen,**  
**0.3ml Hyaluronic Acid Pen**



### More Images



## Product Description

### Description:

1. Please Note: This 2 in 1 Hyaluron Pen is only a tool you need when you work. It does not contain Hyaluron acid. You need to buy hyaluron acid elsewhere. To avoid wasting expensive products, we recommend that you learn hyaluron pen training and practice before using it. Before using it, please practice according to the operation steps in the manual.
2. Painless: no needle injection design, avoid contact with subcutaneous nerve tissue, can greatly reduce the pain of injection.
3. Fast: Fast injection in less than 1/3 of a second, micro orifice, and the automated spring pressure make injections virtually painless, whether it is suitable for home use or beauty salons.
4. Function: Deeply moisturizes and moisturizes the skin, allowing the skin to quickly fill, improving wrinkles and dull skin.
5. Atomization type: The high pressure is generated by the pressure device, and the liquid is pushed out from the micro-hole at the top of the front end to form a 0.17mm ultra-fine liquid column, which instantly penetrates the skin and reaches the skin.

### What is the hyaluron pen?

Water mesotherapy is to inject hyaluronic acid into the skin. Hyaluronic acid contains moisture, 200-300 times as rich as one's body. So that it makes the skin smooth, soft, bright, and hydrated after the injection. To keep the skin hydrated, ladies always choose varied skin care products. Skincare products only hydrate moisture on the cutin. So smearing those products cannot deeply provide moisture to the skin. As we know, acne, wrinkles, and inelasticity appear if the skin lacks moisture. Water mesotherapy will solve all the problems because of lack of moisture, such as sallow skin, color spots, and wrinkles... And it will supply moisture to the dermis. Therefore water mesotherapy will work well on whitening, moisturizing, and renewing inelasticity.

### Principle of operation:

Spring and piston, which allow, under pressure, quickly and evenly inject drugs under the skin through a hole with a diameter of 0.1 mm, which is two times smaller than the needle thickness for injection!

With the help of the energy of mechanical action, the spring and the special piston in the apparatus allow, under pressure, quickly but gently pushing the drug stream through a microscopic hole into the subcutaneous tissue.

Drugs are absorbed much better and faster than with the introduction of a needle, which was confirmed in a series of studies and tests conducted in different countries of the world.

The ideal solution for those prone to scarring and subcutaneous clots.

### Advantages of Hyaluron Pen:

Painless injection

Injection without a needle eliminates the risk of infection from the needle

Does not leave scars and puncture marks on the skin tissue with a needle

Better and faster absorption of drugs

Strictly subcutaneous injection – eliminating the risk of penetration of the drug into the muscle

Fast, easy, and safe

Separate and economical dosing of the drug



### What is the hyaluron pen?

Water mesotherapy is to inject hyaluronic acid into skin. Hyaluronic acid contains moisture, 200-300 times as rich as one's body. So that it makes the skin smooth, soft, bright and hydrated after the injection. In order to keep the skin hydrated, ladies always choose varied skin care

products. Actually, skin care products only hydrate moisture on the cutin. So that over smearing those products cannot deeply provide moisture for the skin. As we know, acne, wrinkle and inelasticity would appear if skin lacks moisture. Water mesotherapy will solves all the problems by reason of lacking moisture, such as sallow skin, color spot, wrinkle... And it will supply moisture for dermis. Therefore water mesotherapy will works well on whitening, moisturizing and renewing inelasticity.



**Principle of operation:**

Spring and piston, which allow, under pressure, quickly and evenly inject drugs under the skin through a hole with a diameter of 0.1 mm, which is two times smaller than the needle thickness for injection!

With the help of the energy of mechanical action, the spring and the special piston in the apparatus allow, under pressure, quickly but gently pushing the drug stream through a microscopic hole into the subcutaneous tissue.

Drugs are absorbed much better and faster than with the introduction of a needle, which was confirmed in a series of studies and tests conducted in different countries of the world.

The ideal solution for those prone to scarring and subcutaneous clots.

# Hyaluron pen accessories



## Advantages of Hyaluron Pen:

- \* Painless injection
- \* Injection without needle eliminates the risk of infection from the needle
- \* Does not leave scars and puncture marks of the skin tissue with a needle
- \* Better and faster absorption of drugs
- \* Strictly subcutaneous injection - eliminating the risk of penetration of the drug into the muscle
- \* Fast, easy and safe
- \* Separate and economical dosing of the drug



Dermal Lip Fillers



Hyaluron Pen

Reasons to **Not Use** the Hyaluron Pen



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA