



shijiazhuang tiantian technology co.,ltd
hadermal-filler.com

HYALURONIC ACID
DERMAL FILLER
INJECTABLE HYALURONIC ACID

PDO Thread Lift

• The thread is made of high quality PDO material, which is safe and effective for lifting the skin.
• It is made of 100% biodegradable material, which is naturally absorbed in the skin and does not cause any side effects.

Thread Hilos Tensores De Pdo Barbed Cog 3D Sharp Needle Pdo Wires Lift

Our Product Introduction

for more products please visit us on hadermal-filler.com


Basic Information

- Place of Origin: china
- Brand Name: DBX
- Certification: CE
- Model Number: PDO PLLA PCL
- Minimum Order Quantity: 10pcs
- Price: 0.5\$-5\$
- Packaging Details: Auminium bag
- Delivery Time: 3-5days
- Payment Terms: T/T, Western Union, MoneyGram, paypal
- Supply Ability: 500000box/month



Product Specification

- Product Name: Thread Facelift
- Materials: Pdo PCL PLLA
- Needle Size: 18g 19g 20g 21g 23g 25g 26g 27g 29g 30g
- Thread Size: 13mm 25mm 38mm 50mm 60mm 70mm 90mm 100mm 120mm
- Thread Type: Mono/ Screw /Tornado/ Multi /Mesh/Nose Thread/Eye Thread
- Thread Type 2: 4D 5D 6D 8D 12D COG/Press Cog /Fishbone Cog Thread
- Packing: Aluminium Bag
- Shelf Life: 2 Years
- Highlight: PDO barbed thread lift, 3D sharp PDO wires, PDO cog thread needles



More Images



Product Description

Thread Hilos Tensores De Pdo Barbed Cog 3D Sharp Needle 21g 60mm Pdo Wires Lift

PDO Thread Lift is the latest and revolutionary treatment for skin tightening and lifting as well as V-shaping the face. These threads are made of PDO (polydioxanone) material which is similar to the threads used in surgical stitches. The threads are absorbable and therefore will be reabsorbed in 4-6 months time leaving nothing behind but the skin structure created which continues to hold for another 15-24 months.

Areas that can be treated include lifting of eye brow, cheeks, corner of the mouth, nasolabial folds fold and neck. With correct placement of the threads, you shall notice a more defined jawlines and the face appear more "V" shaped. Since absorbable sutures are used, there will not be any foreign body in the skin after 6 months.

After cleansing and sterilization of the face, anesthetic in the form of cream or direct injection can be given to reduce discomfort sensation. The doctor will choose the most appropriate type of threads and place it accordingly at various part of your face. The procedure takes about 30 minutes.

Product name	facelift thread
Materials	PDO PCL PLLA
Items	Single, Double, Screw, Knot, COG Thread , Blunt & sharp needle kit
Thread Size	PDO thread USP 5/0,6/0,7/0
Thread Length	30MM-190MM
Needle Length	25MM-150MM
Mono	widely used in skin recovery and local compact function for various body and face areas
Mono screw	Single screw is a single helix, and the pdo screw is raised to the ideal part, which promotes the skin's tight and natural anti-wrinkle face. facial
Tornado	moisturizing, wrinkle, whitening, tender skin, shrink pores, and enhance the profile
Tornado Screw	with 2 strands of screw shaped sutures can hold skin more tightly and stimulate collagen twice more than normal Screw type.
cog	COG with 360 ° directional cogs acts to hold the skin like hook when the thread is pulled. It can be applied in the area of face for lifting.

They're made from a colorless polyester that breaks down in your body after about 6 months.

The presence of these sutures in your skin triggers cells in your body called fibroblasts to produce more collagen .Collagen is the protein that gives your skin its structure and elasticity. Loss of collagen is one of the main causes of aging skin.

PDO threads can be further divided into three categories:

PDO mono threads. Smooth sutures that help rejuvenate your skin by stimulating collagen production.








PDO cog threads. These threads have barbs that latch into your skin like small fishhooks, to provide support and lift parts of your face.




PDO screw threads. Made up of one or two intertwined threads, these are used to help restore volume to sunken parts of your skin.


Other types of threads

PLA and PCA threads are newer than PDO. They last longer in your body and tend to stimulate more collagen production. PLA threads take about 12 months to be absorbed and PCA takes about 12 to 15 months.

Each type of thread is best suited for a particular function. PDO threads are better at while PLA and PCA threads are better at lifting sagging parts of your skin.

	Specifications 26G * 30/38/60/90mm 27G * 25/38/50/60mm 29G * 25/38/50/60mm 30G * 25/38mm		Specifications 26G * 25/38/60/90mm 27G * 25/38/50/60mm
	Specifications 26G * 30/38/60/90mm 27G * 25/38/50/60mm		Specifications 26G * 30/38/60/90mm 27G * 25/38/50/60mm 29G * 25/38/50mm 30G * 25/38mm
	Specifications 26G * 30/38/60/90mm 27G * 25/38/50/60mm		Specifications 18G * 60mm 19G * 38/50/60mm 21G * 38/50/60mm
	Specifications 26G * 30/38/60/90mm 27G * 25/38/50mm		Specifications 21G * 38/50/60mm 23G * 38/50mm
	Specifications 18G * 100/120mm 19G * 100/120mm		Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm
	Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm		Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm
	Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm		Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm
	Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm		Specifications 18G * 60/100/120mm 19G * 38/60/100/120mm 21G * 38/60/90/100mm 23G * 38/60/90mm

	Specifications Needle 100mm 20G * Thread 210 - 560 mm		Specifications 30G * 25/38mm
Press Cog (Whale teeth shape) 	Specifications 18G * 100/120mm 19G * 100/120mm	Press Cog (Fish bone shape) 	Specifications 18G * 100/120mm
Mesh L Needle 	Specifications 18G * 60mm 19G * 38/50/60mm	Multi Thread L 	Specifications 19G * 38/50/60mm 21G * 38/50/60mm 23G * 38/50/60mm



It is premium lifting thread made of PDO which is absorbed after proliferating collagen by stimulating fibroblast. Any of you can experience overall face contouring, skin elasticity. It brings natural and dramatic effect to face and body after simple and safe treatment. The superior quality of the Ultra-thin wall Needle allows the treatment easy, atraumatic and painless.



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA