

Laboratory Equipment 6PCS 20ml/10ml Prp Centrifuge Tubes Centrifugal Separator Blood Plasma Centrifuge

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: Laboratory Equipment
- Model Number: Laboratory Equipment
- Minimum Order Quantity: 1 piece
- Price: \$80-\$100
- Packaging Details: Box
- Delivery Time: 1-3 Days
- Payment Terms: T/T, Western Union, Paypal, Card, Bitcoin
- Supply Ability: 10000000 boxes /day



Product Specification

- Product: Laboratory Equipment 6PCS 20ml/10ml Prp Centrifuge Tubes Centrifugal Separator Blood Plasma Centrifuge
- Type: Laboratory Equipment Prp Centrifuge Tubes
- MOQ: 1 Piece
- Highlight: 20ml PRP centrifuge tubes,
10ml blood plasma centrifuge,
laboratory centrifugal separator tubes



More Images



for more products please visit us on hadermal-filler.com

Product Description

Laboratory Equipment 6PCS 20ml/10ml Prp Centrifuge Tubes Centrifugal Separator Blood Plasma Centrifuge



Steps for usage:

1. Put the machine on a level and strong table and keep the working room clear, dry and ventilated.
2. Indoor temperature 5-32°C is ok. Open the power supply and make sure the indicator lamp is light.
3. Open the cover of the machine, clear the rotor and container.
4. Put the test tube (with same weight) into the machine.
5. Put down the cover of the machine, adjust the timer as your requirement.
6. Switch ON the power. Adjust the governor as you requirement.
7. Before you are going to stop the work, you should adjust the governor to "0" first, and then clean the container and rotor.



Laboratory centrifuge

800 Series electric centrifuge is widely used in hospitals, chemical and Biological Chemistry Laboratory for serum, plasma, immune qualitative analysis experiment. Has beautiful shape, large capacity, small volume, complete function. Its performance is stable, the speed is adjustable and can automatically adjust the balance, low temperature, high efficiency and wide applicability and advantages.

Technical Specification:

Model No.:	800	800B/D
Max Speed:	4000 rpm	4000 RPM
Capacity:	20(ml)×6	20(ml)×6
Max RCF:	1790xg	1790xg
Power Source:	AC220V 50HZ/AC 110V 60HZ	AC220V 50HZ/AC 110V 60HZ
Timer:		0-60 Min
Packing Size:	29*23*26cm	29*23*26cm

Instruction:

1. Please read the operational manual carefully before using the machine.
2. Place empty centrifugal tube evenly before working, after working 1-2 minutes at the maximum speed, then begin to work if there's no problem.
3. Do not touch the rotor while it is running.
4. Do not operate the machine if the rotor is not stably loaded.
5. Please ask the technician to solve the problem if the wire and the plug are damaged or the electric appliances are in trouble. It will cause quality problem if not using recommended appliances by the manufactured factory.
6. Do not use the machine directly in the sunlight.
7. The value of rated parameters of the power plug offered by the consumer should be more than that of the machine, and it should have well earth connection.
8. Keep it away from high temperature machines.
9. Prevent it from strong collision.
10. Pull out the plug from the socket before changing the fuse.
11. Do not open the centrifugal cover of the machine in order to prevent any damage by outer sundries.
12. Keep it out of reach of children.
13. Do not clean the machine with metal material instead of soft cloths with alkali solution.



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA