



Revolax Hyaluronic Acid Filler Korea With Lido Deep 1.1ml

Basic Information

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity: 1box
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms: T/T Credit Card , Western Union,
- Supply Ability:



Product Specification

| Highlight: | Revolax Hyaluronic Acid Filler, 1.1ml Hyaluronic Acid Filler, |
|--------------------------------------|--|
| Application: | Lips /Nose /Wrinkle |
| Packing Include: | 2 X 27G Needle |
| Lido Concent: | 3mg/ml Lido |
| Concentration: | 24 Mg/ml Cross-linked Hyaluronic Acid |
| • Type 3: | Revolax Sub-Q |
| • Type 2: | Revolax Deep |
| • Type 1: | Revolax Fine |
| • Product Name: | Revolax |
| | |

Korea

CE

35

box

5-8days

MoneyGram

5000box/month

Revolax

Fine/Deep/Sub-Q



Our Product Introduction

More Images





Lidocain Deep Hyaluronic Acid Filler

Product Description

Revolax Deep Fine Sub-Q (1 x 1.1ml) Hyaluronic Acid Dermal Filler

Description :

Revolax Sub-Q offer a smooth and pain-free experience in augmenting the facial features of the patient. The gel contains 0.1ml of Lido , which means any immediate discomfort, typical for a dermal filler injection procedure is not present.

The Revolax product is made from purified Hyaluronic Acid and its gel is monophasic in structure – an advantage, when compared to most filler brands, as its structure is more durable to the enzymes of the body and thus, its effects are more lasting. Revolax Sub-Q with Lido will last 6 months more, on average, compared to most other sculpting filler brands.

The filler has a very high concentration of Hyaluronic Acid – 24mg/ml. Considering how powerful it's effects on the skin are, it is no wonder this is an excellent anti-age product. While using synthetic HA, the substance itself is natural to the human body, with capability of retaining many times it's own weight of water, is completely safe and biodegradeable.

What areas is Revolax Sub Q best injected in?

The primary focus of thick fillers is enhancing the shape of the various facial features of the patient. Revolax Sub-Q with Lido is no exception and it can help doctors improve the patient's:

Jawline Cheek volume Cheekbones Chin Nose (rhinoplasty)

Intended for:

Deep wrinkles Nasolabial folds Augmentation of the cheeks, chin and lips

Specification :

| Product name | Revolax |
|----------------|--------------------|
| Туре | Fine /Deep/Sub-Q |
| Concent | 24mg/ML |
| Lido concent | 3mg/ML |
| pack consists | 1 x 1.1ml syringe |
| | 2 x 27G needle |
| Delivery time | within 4-6days |
| Delivery Terms | DHL,FEDEX,UPS |
| | Jawline |
| | Cheek volume |
| Application | Cheekbones |
| | Chin |
| | Nose (rhinoplasty) |

| REVOLAX" FINE REVOLAX" FINE with Lideceline | Fine Lines Crow's feet Glabella Lines Forehead Lines Earlobe Perioral wrinkles | |
|---|--|---|
| REVOLAX [®] DEEP REVOLAX [®] DEEP with Addressive | Deep Lines / Augmentation Nasolabial Lines Glabella Lines Forehead Lines Laugh Lines Cheek/Chin/Forehead augmentation Lip contouring Nose Augmentation | |
| REVOLAX" SUB-Q REVOLAX" SUB-Q with lideoing | Very Deep Lines / Augmentation Nasolabial Lines Laugh Lines Cheek/Chin/Forehead augmentation Lip contouring Nose Augmentation Breast Augmentation | 9 |

