



Factory Directly Sale Skin Care Beauty Double-Headed Hyaluronic Acid Pen

Our Product Introduction

Basic Information

- Place of Origin: china
- Brand Name: hyaluron pen
- Certification: ce
- Model Number: 0.3ml/0.5ml
- Minimum Order Quantity: 1
- Price: usd 30
- Packaging Details: 1 0.3ml ampoule, 1 0.5ml ampoule, 1 pen, 1 box
- Delivery Time: 3-5day
- Payment Terms: Western Union, MoneyGram, T/T, paypal, credit card
- Supply Ability: 100000 pcs



Product Specification

- Packing Size: 22.0 * 12.0 * 7.0 (Cm)
- Applications: Wrinkle Remover And Lip Lifting
- Content: 1 0.3ml Ampoule, 1 0.5ml Ampoule, 1 pen, 1 Box
- Technology: Needle Free High Pressure Atomization



More Images



Product Description

0.3 0.5ml Hyaluron Pen wrinkle removal Hyaluronic Acid injection Lip Dermal filler Pen



Product Description:

Type: Hyaluronic Acid Pen

Feature: Blood Vessels Removal, Skin Rejuvenation, Dark Circles, Wrinkle Remover, Weight Loss, Moisturizer, Whitening, Moisturizer

Certification: Medical CE / ISO13485 / RoHS

material: 304 Stainless Steel

Injection volume: 0.01ml~0.12ml (The min scale 1 frame is 0.01ml, if you want to play 0.05, then down 5 frames, a maximum of 12 frames is 0.12ml the biggest scale is 0.5ml and another style is 0.3ml)

Operation Mode: High Pressure Atomization

DUAL Hyaluron Pen MAX

0.3ml / 0.5ml
Ampoule

3 Pressure
Levels

Scale

0.3ml / 0.5ml
Ampoule Connector

Injection Button

What is the hyaluron pen?

Water mesotherapy is to inject hyaluronic acid into skin. Hyaluronic acid contains moisture, 200-300 times as rich as one's body. So that it makes the skin smooth, soft, bright and hydrated after the injection. In order to keep the skin hydrated, ladies always choose varied skin care products. Actually, skin care products only hydrate moisture on the cutin. So that over smearing those products cannot deeply provide moisture for the skin. As we know, acne, wrinkle and inelasticity would appear if skin lacks moisture. Water mesotherapy will solves all the problems by reason of lacking moisture, such as sallow skin, color spot, wrinkle... And it will supply moisture for dermis. Therefore water mesotherapy will works well on whitening, moisturizing and renewing inelasticity.



Principle of operation:

Spring and piston, which allow, under pressure, quickly and evenly inject drugs under the skin through a hole with a diameter of 0.1 mm, which is two times smaller than the needle thickness for injection!

With the help of the energy of mechanical action, the spring and the special piston in the apparatus allow, under pressure, quickly but gently pushing the drug stream through a microscopic hole into the subcutaneous tissue.

Drugs are absorbed much better and faster than with the introduction of a needle, which was confirmed in a series of studies and tests conducted in different countries of the world.

The ideal solution for those prone to scarring and subcutaneous clots.

Hyaluron pen accessories

Connector head	 0.5ml	 0.3ml
Ampoule head	 0.5ml	 0.3ml
Needle (don't include needle)	 0.5ml	 0.3ml

Advantages of Hyaluron Pen:

- * Painless injection
- * Injection without needle eliminates the risk of infection from the needle
- * Does not leave scars and puncture marks of the skin tissue with a needle
- * Better and faster absorption of drugs
- * Strictly subcutaneous injection - eliminating the risk of penetration of the drug into the muscle
- * Fast, easy and safe
- * Separate and economical dosing of the drug



Dermal Lip Fillers



Hyaluron Pen

Reasons to **Not Use** the Hyaluron Pen



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA