

china

Ce

35\$

box

3days

MoneyGram

50000box/month

Replengen

volume ,deep

T/T, L/C, D/A, D/P, Western Union,

Replengen Volume 1.1MI Hyaluronic Acid Dermal Filler Needle Size 27G

Basic Information

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity: 1box
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms:
- Supply Ability:



Product Specification

• Product Name:	Replengen Deep
 Materials: 	Hyaluronic Acid
• Size:	1ml
• Grade:	Volume
• Grade 1:	Deep
Composition:	HA 24 Mg/ml
• Lido:	0.3%
Needle Size:	27G
 Highlight: 	1.1MI Hyaluronic Acid Dermal Filler, 27G Needle Hyaluronic Acid Dermal Filler



More Images







Product Description

Replengen Volume - Reliable Dermal Fillers 1.1 ml Hyaluronic acid Dermal Filler Product Descrption :

Replengen is a dermal HA-based filler for correction of medium and severe wrinkles. Replengen Deep is used for injection into deep subcutaneous layers. Special technologies are used in the production of Replengen Deep: "Glenzen's GI (Grinding Irregular)" ensures a homogeneous gel structure, while HAMC (Hyaluronic Acid Moisture Correction) technology enables obtaining highly cross-linked HA. The dermal filler has lower rates of decomposition when injected in the epidermis thus increasing the effect duration.

Product name	Replengen Filler
Туре	Volume /Deep
composition	HA 24 mg/ml, Lido 0.3%
Effect	12 to 18 months.
Packing Include	1 syringe × 1.1 ml; 1 needle, 1 cannula per pack

Strengths of Replengen:

the effect is noticeable right after the procedure

improves skin elasticity and moisturizes it

the presence of lido in the composition and homogeneous gel structure make the procedure painless when injecting into deep skin layers the product is made using highest-quality HA

clinically proven efficacy and safety

Replengen Deep is perfect for filling of moderate and deep facial wrinkles. Product is an analogue of UltraFill Deep.

Scope of Replengen Deep:

correction of moderate-to-deep wrinkles (nasolabial fold, forehead and frown wrinkles, lip lines, marionette lines) lip augmentation facial volume augmentation chin correction The aesthetic effect lasts from 12 to 18 months.

Product composition: HA 24 mg/ml, Lido 0.3%

Replengen Deep 1 syringe × 1.1 ml; 1 needle, 1 cannula per pack Needle size: 27G Manufacturer: Reanzen Co., Ltd., South Korea

The product has contraindications. We are not responsible for any negative effects that you may encounter after administering products by yourself. Consulting a doctor or medical professional is recommended prior to any procedure.

Replengen Volume lido 1.1ml

Over time, the skin loses its elasticity, which is visible, among others, in the form of wrinkles of various depths. In order to fill deep wrinkles and prevent new ones, it is worth reaching for Replengen Volume filler.

What is Replengen Volume lido ?

Replengen Volume is an innovative Korean filler based on cross-linked hyaluronic acid. Manufactured by Reanzen Co., Ltd., it is the third in the Replagen series of fillers. It is intended for injection into the middle and deep areas of the dermis. It is produced using two specially developed technologies:

Glenzen's GI (Grinding Irregular) is responsible for the homogeneous, appropriately absorbable gel structure, HAMC (Hyaluronic Acid Moisture Correction) technology is responsible for obtaining highly cross-linked HA.

Thanks to this, the filler decomposes more slowly after injection, which extends the duration of action and prolongs the effect of the treatment.

Application

Thanks to the presence of lido, treatments with Replengen Volume are painless. The use of the product is recommended for the loss of facial skin volume, for filling deep wrinkles, sunken cheeks and drooping skin. It is also great for filling and correcting the lip contour. As a result, the skin not only gains a younger appearance, but also becomes more moisturized, and its condition improves significantly.

The effects of using Replengen Volume

improvement of the face contour, natural rejuvenation effect, reduction of deep wrinkles, delays the formation of new wrinkles and furrows, long-lasting effect of the treatment.

