

# Rejeunesse Deep Hyaluronic Acid Dermal Filler With Lodocain Korea Cross Linked

## **Basic Information**

- Brand Name:
- Model Number:
- Minimum Order Quantity:
- Packaging Details:
- Delivery Time: 4-7
- Payment Terms:
- Supply Ability:
- 1box box 4-7work days L/C, D/A, D/P, T/T, Western Union 1456+box+week

Rejeunesse

1.1ml



## **Product Specification**

- Product Name:
- Function:

• Form:

• Usage:

• Feature:

• Highlight:

- Hyaluronic Acid Dermal Filler Remove Wrinkle Gel Skin Filler/Face Filler,Facial Skineance Wrinkle Filler Cross Linked
- Rejeunesse Deep Hyaluronic Acid Dermal Filler, Hyaluronic Acid Dermal Filler Cross Linked



Our Product Introduction

## More Images









## **Product Description**

## Rejeunesse Deep Hyaluronic Acid Dermal Filler With Lodocain Korea

Non-animal hyaluronic acid & no residual crosslinker.



Rejeunesse is non-animal based, cross linked dermal filler using BDDE for linking agent. Rejeuness product is obtained through our unique BDDE removal process, making it safe and non carcinogenic for humans. Rejeuness contains less than 1.2EU/ml of endotoxins with no BDDE residue detected and it exceeds the highest level of HA purity. Rejeunesse is completely free of all animal products thus reducing the risk of an immunogenic reaction.

-LONG-LASTING.

Immediately visible results after application.

Once injected, HA is gradually broken and absorbed into th body, which is how strongly HA is crosslinked is the crucial to the longevity of the dermal fillers. With unique cross-linking technology, Rejeunesse has longevity over 12months.

### -EFFICACY.

More stable forms of HA, in paticular by 100% croos-linked HA molecoles.

Rejeunesse offers visible results immediately and stable volume effect with minimized

swelling after application.

## -TREATMENT WITHOUT PAIN

Rejeunesse has optimized the product to minimize the pain during and/or post injection for the patients. Rejeunesse enables smoothest and stable injection during procedure by maintaining consistent injection force. Especially, Rejeunesse Lidocain series which combines 0.3% of idocaine(Local anesthetic)minimizes injection pain during the procedure. Thus, it allows patients to experience the most comfortable, almost painless, procedure of dermal fillers.

	with idocaine ( 24mg/ml, 1.1ml )
Rejeunesse FINE	This is used in the correction of fine to medium size depression of the face via mid to superficial dermis injection.
Rejeunesse DEEP	with idocaine (24mg/ml, 1.1ml) This is used in the correction of medium to deep sized depression of the face via mid to deep dermis injection
Rejeunesse SHAPE	with idocain ( 24mg/ml, 1.1ml ) This is used in the correction of deep sized depression of the face via deep dermis to suncutaneous injection

### LONG-LASTING.

Immediately visible results after application.

Once injected, HA is gradually broken and absorbed into th body, which is how strongly HA is crosslinked is the crucial to the longevity of the dermal fillers. With unique cross-linking technology, Rejeunesse has longevity over 12months.

### EFFICACY

More stable forms of HA, in paticular by 100% croos-linked HA molecoles. Rejeunesse offers visible results immediately and stable volume effect with minimized swelling after application.

### **Rejeunesse FINE**

with idocaine ( 24mg/ml, 1.1ml )

This is used in the correction of fine to medium size depression of the face via mid to

superficial dermis injection.

#### **Rejeunesse DEEP**

with idocaine ( 24mg/ml, 1.1ml )

This is used in the correction of medium to deep sized depression of the face via mid to

deep dermis injection

#### - Rejeunesse SHAPE

with idocain ( 24mg/ml, 1.1ml )

This is used in the correction of deep sized depression of the face via deep dermis to

suncutaneous injection

