



# Revolax Fine Deep Sub-Q Hyaluronic Acid Dermal Filler

# **Basic Information**

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity: 1ml
- Price: usd 28
- Packaging Details:

• Supply Ability:

- Delivery Time: 3-5day
- Payment Terms: T/T, Western Union, MoneyGram, credit card

usa

ce

1ml

revolax

one box 1 pcs

,paypal 100000pcs/day



# Product Specification

Product Name:	Hyaluronic Acid Gel,Hyaluronic Acid Dermal Filler,hyaluronic Acid Dermal Injection Filler 10 MI For Skin Tightening,Lip Dermal Filler HA Injection Cosfiller Brand Hyaluronic Acid Gel
• Function:	Remove Wrinkle,anti-aging,breast Enhancement
• Material:	Cross-linked HA,Hyaluronic Acid Filler,Non Animal Cross-linked Hyaluronic Acid
Application:	Breast And Buttock Enlargement,Body Filler,Lip Enhancement,Anti-Wrinkle
<ul> <li>Product Name:</li> </ul>	Revolax
<ul> <li>Packing:</li> </ul>	1*1ml
• Type1:	Fine
• Type2:	Deep
• Type 3:	Sub Q
<ul> <li>Highlight:</li> </ul>	Hyaluronic Acid Dermal Filler,

Cub O Hughurania Aaid Darr

## WhatsApp: 008618132165209 Email: nancy@h&kouxuan.com



# More Images

WhatsApp: 008618132165209 Email: nancy@hbkouxuan.com



### CE Revolax Deep Sub-Q Fine With Lido Derma Filler Injection Hyaluronic Acid

## Revolax Deep/Sub-Q/Fine

### **Discover Our Revolax Dermal Fillers Range**

At Filladerm, you'll find all three types of fillers offered by Revolax. All the fillers are designed to beautify different facial areas. While cross-linked hyaluronic acid is the main ingredient of all of them, the main difference is in the viscosity. Another thing to mention is these soft fillers are available with

**Revolax Fine** 

Revolax Fine is the thinnest dermal filler in its range. Its complex viscosity is 80 MPa which makes it the perfect choice for treating superficial lines, crow's feet, glabellar lines, neck wrinkles, and lip augmentation. The results of Revolax Fine last for 8 - 12 months. Revolax Deep

Revolax Deep is also based on pure hyaluronic acid with a concentration of 20 mg/ml. Revolax Deep is thicker than Revolax Fine. Its complex viscosity is 160 MPa.

The effect of Revolax Deep lasts for 12 - 18 months. It's used in the treatment of nasolabial folds and augmentation of the cheeks, chin, and lips. Revolax Sub-Q

Revolax Sub-Q is the thickest filler within its range in terms of the viscosity of the hyaluronic acid gel. Its complex viscosity is 280 MPa. This thickness of the gel makes it a good choice when it comes to treating the augmentation of the cheeks, chin, nose, and jawline. The result should last 12 - 18 months.

#### Where to use Revolax HA Gel?

Revolax Deep can be used by the lips, cheeks, face contouring, and marionette lines. Revolax Fine and the Revolax Deep can be used on the Nasolabial folds, deep wrinkles, Nose, Cheeks, and chin.



#### spec

Desident		
Product	Hyaluronic Acid Gel, Hyaluronic Acid Dermal Filler, hyaluronic acid dermal injection filler 10 ml for skin tightening, Filler for skin, Lip	
name	dermal filler HA injection Cosfiller Brand Hyaluronic acid Gel	
Function	Cross Linked, safety, easy to use	
Cross-linked		
agent	DUDE	
Packing	Blister Package, 1syringe/box, 1pc/box, One Prefilled Syringe with Two Needles.	
Material I	Pure Hyaluronic Acid, Cross-linked Sodium Hyaluronate, Cross-linked HA, Hyaluronic acid filler, Non Animal Cross-linked	
	Hyaluronic Acid	
Application	breast and buttock enlargement, Body Filler, Lip enhancement, Neck, Anti-Wrinkle	
Туре	fine,derm.deep.ultra deep	
Ha		
concentratio 24mg/ml, 20mg30mg/ml		
n		
Particle size	0.28~0.5mm, 0.1mm-2.0mm	
Duration	6~24 months, 6-18 months, 6-18months, 8~12 Months	

Revolax Fine: This is the Hyaluronic Acid with the lowest viscosity, it is the best for smoothing small lines on the face, neck and body. It eliminates the first signs of ageing and prevents wrinkles deepening. Works like hydration for the skin. It lasts between 6 to 9 months.

Revolax Deep: is a hyaluronic acid gel with medium viscosity. It efficiently fills in wrinkles and folds on the skin and works well on the deep dermal layers

Duration period of 8 to 10 months.

Revolax SUB Q: is designed to fill the volume of areas on the face. It is the best for deep wrinkles, skin folds, huge lips with a lot of volume and facial contouring. This hyaluronic acid gel is the thickest, compared to all the other Revolax products



#### How does the treatment work?

The treatment can be done in 15 minutes but it can also take longer, depending on the area to be treated. The professional in charge of your treatment injects the product using an ultra-thin needle under the surface of the skin. For fine lines, the product is injected at the surface of the skin. For marked wrinkles, the injection is made at a deeper level. Please note that the injection of Juvéderm Ultra requires a medical gesture that can only be carried out by a doctor (doctor specialised in dermoesthetics, dermatologist or plastic surgeon).



