

Disposable Micro Cannula Needle Blunt Tip Micro Blunt Cannula for Hyaluronic dermal filler

Basic Information

- Brand Name:
- CANNULA NEEDLE
- Certification:
- Minimum Order Quantity: 10 pcs
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms:
- Supply Ability:
- 0.8 \$ Sterile packaging 1 to 5 days T/T, Western Union, Credit Card
- y: 50000 pcs/month

CE



Product Specification

- Product Name:
- Usage:
- Feature:
- Packing:
- Size:
- Material:
- Blunt Tip Micro Blunt Cannula Hyaluronic Acid Dermal Filler Injection Safety Sterile Packaging 14g To 34g, 25mm To 70mm
- Sus304



More Images




Product Description

Factory Price Top Quality 304 Stainless Steel Blunt Cannula All Sizes 18g-30g Blunt Needle



A blunt cannula allows blood vessels to be 'pushed aside' as the microcannula moves around the tissue under the skin instead of piercing the blood vessels as with traditional needles. Safety from injuring nerves Safety from injuring blood vessels. Safety from injuring glands and ducts. Better Results: you can get a smoother and wider fanning result. Better Results: you can usually have less bruising with your patients. Better Pain management: Patients usually feel less pain with this technique.





Flexible Cannula FLEXIBLE TO STIFF TO SUIT ALL NEEDS Flexible cannulas for superficial & stiffer cannulas for deep injections. Blunt tip DOES NOT CUT THROUGH NERVES OR VESSELS Dermasculpt's blunt tips slide under the skin without causing any trauma. Painless

DermaSculpt decreases bruising. Safer than needles, it allows for less traumatic injections.

