


## Beauty Injection Micro Cannula Blunt Tip Needle 14g 15g 16g 18g 22g 38mm 50mm 60mm 70mm

Our Product Introduction

for more products please visit us on hadermal-filler.com

### Basic Information

- Place of Origin: China
- Brand Name: Micro Cannula Needle
- Certification: CE
- Model Number: 18G 21G 23G 25G 27G 30G
- Minimum Order Quantity: 10
- Price: 0.7
- Packaging Details: Sterile package
- Delivery Time: 1-3Days
- Payment Terms: T/T, Western Union
- Supply Ability: 80000




### Product Specification

- Product Name: Micro Cannula Needle
- Usage: Single-Use
- Type: Blunt Tip
- Size: 25mm 38mm 50mm 70mm
- Disinfecting Type: Far Infrared
- Material: Stainless Steel
- Function: Injection & Puncture Instrument
- Warranty: 2 Years
- Spe: 14g 16g 17g 18g 19g



### More Images



### Product Description

#### Disposable Stainless Blunt Tip Needle Micro Cannula Needle 18G 21G 27G

##### What Is A Blunt Cannula?

Cannulae are blunt, which mean that they are less likely to pierce a through a vein or artery. This results in less swelling, bruising and intravascular injection of filler. The practitioner has more feel about what structures are being passed by the tip of the cannula.

##### What Is A Blunt Plastic Cannula Used For?

The BD blunt plastic cannula is a truly universal cannula designed to be used with all types of split septum IV injection sites and vials for needle-free access.

##### How Do You Use A Blunt Tip Cannula?

A needle is required to make an entry point for the cannula that just passes through the skin. The cannula is then used to pass through the entry point to the areas beneath the skin. The blunt-tipped cannula can re-enter the same entry point again without having to hurt the patient with another needle puncture.

|       |                      |
|-------|----------------------|
| Name  | blunt cannula needle |
| size  | 15g 16g 18g 19g      |
| size2 | 38mm 50mm 60mm 90mm  |


Our Product Introduction



## Micro-Cannula

Micro-cannula is an improved version of the traditional hypodermic needles used to administer injection treatments.

### Product specifications




| Site | Length(mm) | Site | Length(mm) |
|------|------------|------|------------|
| 18G  | 50         | 25G  | 38         |
|      | 70         |      | 50         |
|      | 100        |      | 70         |
| 21G  | 50         | 27G  | 38         |
|      | 70         |      | 50         |
| 22G  | 50         |      | 70         |
|      | 70         | 38   |            |
| 23G  | 38         | 30G  | 40         |
|      | 50         |      | 50         |
|      | 70         |      | 70         |

## Product features

**A/** Flexible cannulas for superficial & stiffer cannulas for deep injections.



Blunt tips slide under the skin without causing any trauma. **/B**





Specifically designed for atraumatic intradermal injections of fluids.



## Product advantages



## Product shows



shijiazhuang tiantian technology co.,ltd

008619931162089

info@tiantianbike.com

hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA