



# Audrey 10ml Injectable Buttock Enlargement Filler Injection M Audrey

## **Basic Information**

Place of Origin: KoreaBrand Name: fillerCertification: CE

Model Number: 1.1ml 10ml
Minimum Order Quantity: 1box
Price: 78
Packaging Details: box
Delivery Time: 1-3days

Payment Terms: T/T, Western Union, MoneyGram,credit card

• Supply Ability: 50000000box/month



## **Product Specification**

Product Name: Filler 10mlType: 10 Ml X1pcs

• Packing: 1 X 1.1 MI Syringe/box

Typically Performed: 15-30 MinsLast Time: 6-18 Months

• Lidoos: 0.3%

Applications: Facial Lips Nose Check

Needles: 27GType 1: LYFTDelivery Time: 1-2days

• Highlight: Audrey Injectable Buttock Enlargement Filler,

M Audrey Injectable Buttock Enlargement Filler, 10ml Injectable Buttock Enlargement Filler





## **Product Description**

#### **M** Audrey


| product    | dermal filler hyaluronic acid lips facial |
|------------|---|
| Туре       | 1.1 10                                    |
| size       | 1.1ml x 1                                 |
| effection  | dermal filler                             |
| how to use | injection                                 |
| usage      | lips check facial nose filler             |
| usages     | dermal filler injection                   |
| MOQ        | 1 piece                                   |
| eccects    | Anti aging wrinkels                       |

AUDREY FILLERS are manufactured by Ostin Bio, Co., Ltd using hyaluronic acid derived from non-animal origins with high purity; hence it is biocompatible when injected into the body. AUDREY FILLERS are composed of highly purified and internationally recognized derivatives of hyaluronic acid that constantly maintain elasticity in the skin.

AUDREY FILLERS are developed with MSFT or Multiple Swelling Formulation Technology, a core method patented by Ostin Bio. It is a method where the hyaluronic acid in the dermal filler uses minimal cross-linking agents. It uses a blend of gels with different swelling ratios to increase hydration and safely maintain volume once injected into the skin. This innovative procedure gives AUDREY FILLERS an excellent volume effect, outstanding tissue adhesion, and high cohesiveness. It also minimizes free radical formation by significantly reducing proteins that may cause hypersensitivity reactions. Upon clinical trials, AUDREY FILLERS shows low unit levels of endotoxin, thus holding an outstanding safety profile.

AUDREY FILLERS come in different variants according to their consistency. The hard variation of the AUDREY Filler is manufactured with a higher viscosity and more volume than the soft type, thus resulting in a more robust shaping capacity after injection. It somehow prevents deformation and can be effective for a longer duration.











#### **Description**

1ml with s an efficiant dermal filler for correcting moderate to severe lines and wrinkles as well as plumping and shaping your patient's lips. Product is a sterile, transparent gel of stabilized hyaluronic acid of non-animal origin with the addition of 0.3% hydrochloride. 1ml with is supplied in a glass syringe with a Luer lock fitting. Only authorized aesthetic professional in accordance with local regulations should administer with injections.

#### The packaging content:

- 1 prefilled syringe 1,0 ml
- 2 x needles
- Tracebility labels

version is suitable for improvement for fine lines, rinkles and folds. It is highly recommended to use for lips enhancment. The brand is made using patented Non-Animal Stabilized Hyaluronic Acid (NASHA) Technology, which helps retain the natural way that hyaluronic acid molecules connect to each other. This results in a product that produces much longer-lasting results. As well, this filler comes with added, an anesthetic for better patient experience during injection.

Benefits:

| Model NO.     | radiesse |
|---------------|----------|
|               |          |
| Delivery Time | 3-5day   |
| Specification | 1ML      |
| Origin        | Korea    |

| Efficacy            | Wrinkle Removal & Skin<br>Beauty |
|---------------------|----------------------------------|
| Transport Package   | Box                              |
| Trademark           | revolax                          |
| Production Capacity | 1000000                          |



injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject Revolax with into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject Revolax with into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.









hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA