



Cross-Linked Elasty for Long-Lasting Hyaluronic Acid Dermal Filler

Our Product Introduction

for more products please visit us on hadermal-filler.com

Basic Information

- Place of Origin: Korea
- Brand Name: Elasty
- Certification: CE
- Model Number: D G F
- Minimum Order Quantity: 1box
- Price: 38USD
- Packaging Details: box
- Delivery Time: 5-8 work days
- Payment Terms: T/T, Western Union, MoneyGram paypal ,credit card
- Supply Ability: 1000box/month



Product Specification

- Product Name: Elasty
- Type 1: ELASTY Plus D
- Type 2: ELASTY Plus F
- Type 3: ELASTY Plus G
- Content: HA 24mg/mL
- Lido: 0.3%
- Packing Size: 1x1ml
- Whatsapp: 008613315976231



More Images



Product Description

Product Description

ELASTY

HA 24mg/mL

Elasty is a dermal filler contains hyaluronic acid (HA) at a concentration of 24mg/ml, with the addition of lidocaine. Elasty is known for its excellent cohesion and high elastic modulus.

ELASTY Plus with Lidocaine hyaluronic acid dermal fillers are top-of-the-line thanks to their patented Preserved Natural Entanglement Technology (PNET), providing enhanced high elastic modulus and excellent cohesion. These fillers boast high elasticity, high cohesion, and high convenience.

With minimal modification, this dermal filler is highly pure and preferred by consumers. Its high viscosity guarantees maximum satisfaction, making it a perfect choice for quality-focused individuals. High viscosity ensures a proper "fill," achieving optimal results for patients. Additionally, the inclusion of lidocaine minimizes discomfort during the procedure.

D plus is a dermal filler with medium viscosity. It efficiently fills in wrinkles and folds on the skin and works well on the deep dermal layers. It has been manufactured by cross-linking 3D technology and its high volumetric properties enable resistance to biodegradation.



Product advantages:

High stage of purification of hyaluronic acid;
Minimal viscosity facilitates an easy injection of the gel;
Cross-linked fibers ensure a long-lasting correction effect;
Ideal for a thin skin and for sensitive skin areas;
Painless procedures.
High stage purification of hyaluronic acid
Minimal viscosity facilitates seamless injections
Cross-linked fibers ensure a long-lasting effect
Ideal for thin skin and for sensitive areas
Lido enables painless procedures

Product purpose

D plus Filler allows filling in moderate to deep wrinkles around eyes, on the forehead, between the eyebrows, cheeks, marionette lines on other areas of face. It gently enlarges natural volume of tissues.

Description and composition

The product is released in 1ml syringes filled with gel containing hyaluronic acid (24 mg/ml)

The Hyaluron Pen is a revolutionary device that uses pressure to inject Hyaluronic Acid into the skin – with no needles, no pain, and no downtime. The Hyaluron Pen is a no-needle, no-pain treatment that creates plump lips, fixes

smile lines, and create instant, long-lasting results. Not only is it needle-free and pain-free, but it also creates beautiful results in just minutes. Made in South Korea.

Application Area:

- Correction of the medium and deep wrinkles
- Correction of the nose and nasolabial folds
- Chin augmentation
- Forehead wrinkles smoothing
- Correction of glabella wrinkles

Main Ingredients:

HA(hyaluronic acid): 24mg/ml

Elasty G – dermal filler based on highly purified hyaluronic acid, developed by PNET (Preserved Natural Entanglement Technology). This technology is unique in increasing the cohesion coefficient of the molecules and, consequently, a stronger structure and longer lasting results. The product is used for the correction of medium and deep wrinkles and folds, volumizing and contouring.

Elasty G does not contain lido

Advantages of Elasty G

- High viscoelasticity
- Pronounced effect after the first treatment.
- Purified HA of non-animal origin
- Minimal risk of side effects, swelling, and allergic reactions
- Cross-linked HA fibres
- Suitable for thin and sensitive skin

Scope of Elasty G

- Correction of the fine lines and wrinkles
- Smoothing of the forehead wrinkles
- Nasal labial folds treatment
- Correction of neck wrinkles and wrinkles of periorbital area
- Lips augmentation











Composition (main ingredients)

HA 24 mg/ml

Duration

The effect of the procedure lasts from 8 to 12 months.

CASE STUDIES — ELASTY GRAND PLUS

 BEFORE	 6mL for cheeks, jawline, chin 6mL = 3 boxes Elasty Grand Note filler along jawline tightens jowls and lifts nasolabial folds	 BEFORE	 2mL for cheeks & chin 2mL = 1 box Elasty Grand
 BEFORE	 6mL for cheeks & jawline 6mL = 3 boxes Elasty Grand note filler creates neck lift	 BEFORE	 6mL Elasty Grand + Revolax Fine 6mL = 3 boxes Elasty Grand Elasty Grand to enhance cheeks, define jawline, restore chin Revolax FINE for under eye circles and bottom lip. Note how filler lifts double-chin and sagging around mouth & nose.
 BEFORE	 1mL for lip augmentation 1/2 syringe (.5mL) used for upper lip and .5mL for lower lip		



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA