shijiazhuang tiantian technology coultd

Liporase Hyaluronidase Is An Enzyme That Dissolves Hyaluronic Acid

Basic Information

Place of Origin: KoreaBrand Name: LiporaseCertification: CE

Model Number: Hyaluronidase

Minimum Order Quantity: 1 box
Price: 50UD
Packaging Details: box
Delivery Time: 5-8days

Payment Terms: T/T, Western Union, MoneyGram

• Supply Ability: 5000box/month



Product Specification

Product Name: Hyaluronidase
Other Name: Liporase
Each Vial Contains: 1,500 I.U
To Be Mix Saline: 3ml

Application: Dissolve Dermal Fillers

Package: 10vials/boxHyaluronidase (EP): 1500 IULactose Hydrate: 13.3 Mg



More Images







10 Vials Liporase Dissolve Hyaluronic Acid Dermal Fillers Hyaluronidase

Specification:

Product name	Liporase Hyaluronidase
Hyaluronidase (EP)	1500iu
Lactose Hydrate	13.3mg
packing	10 vials/box
Application	dissolve the hyaluronic acid
Usage	adding Hyaluronidase to lipodissolving injection
1 Vials	Mix 3cc of saline
1 vials	disslove 3cc of hyaluronic acid

Description

What is Liporase Hyaluronidase and what does it do?

This product contains a specially designed protein enzyme that quickly and effectively dissolves Hyaluronic Acid (HA) and other fluids by speeding up the rate at which they are absorbed by the body. Liporase Hyaluronidase is a risk-free product that rarely leads to adverse effects most of which are a result of injection site reactions and not the product itself. The main use of this tissue modifier are to dissolve HA dermal fillers and similar Hyaluronic acid based gels.

Simply & swiftly remove fillers:

Liporase starts breaking down the Hyaluronic Acid as soon as it is injected into the target area. Depending on the HA concentration, effects become visible after 24~48 hours.

Boosts the effects of fat dissolving injections:

Helps the fat destroyed by fat dissolving injections to be absorbed by the body, boosting the effects.

Hydrolysis of HA fillers

Liporase Hyaluronidase swiftly breaks down Hyaluronic Acid of a synthetic origin. This means that it can successfully reverse the effects of a filler injection gel without causing your body to lose its naturally-produced HA. This product is suitable for the following cases:

Uneven distribution of the HA gel

Overfilling

Incorrect placement of the dermal filler

Accidental intraarterial injection

To dissolve 3cc of Hyaluronic Acid, mix 1 vial of this solution with 3cc of saline water and inject it directly in the treated area

Process of adding Hyaluronidase to lipodissolving injection:

Mix 5cc of saline to 1 vial of Hyaluronidase, it could be add to a 10cc lipodissolving Vline injection. If it is 30cc, add 2 vials of Hyaluronidase. Hyaluronidase is mixed with Vline injection, it is injected along with lido into the fat layer evenly.

Massage the area where it was treated for 15 minutes each day and 7 days in a roll, apply hot packing after massage.

Caution:

There have been a few reports that some people may have an allergy to hyaluronidase, which happens to those who are allergic to bee stings.





shijiazhuang tiantian technology co.,ltd

008619931162089





hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA