# Eptq S100 S300 S500 Hyaluronic Acid Injection Hyaluronic Acid Dermal Filler For Face

## Basic Information

• Place of Origin: Korea Brand Name: e.p.t.q CE · Certification:

s500 s100 s300 Model Number:

• Minimum Order Quantity: 1box 38USD • Packaging Details: box

• Delivery Time: 5-8 work days

• Payment Terms: T/T, Western Union, MoneyGram paypal ,credit card

Supply Ability: 1000box/month



## **Product Specification**

• Product Name: E.p.t.q Type 1: E.p.t.q S100 • Type 2: E.p.t.q S300 Type 3: E.p.t.q S500 • Content: HA 24mg/mL 0.3% • Lido: • Packing Size: 1x1ml

Whatsapp: 008613333015239

Highlight: Face Hyaluronic Acid Dermal Filler,

Hyaluronic Acid Injection Dermal Filler, Eptq S100 Hyaluronic Acid Dermal Filler



## More Images







## **Product Description**

# for more products please visit us on hadermal-filler.com

Our Product Introduction

# Our Product Introduct



Why e.p.t.q. dermal filler?

Dermal fillers help put a little bit of magic into your face. Whether it's a lip plump, smoother skin or youthfulness, dermal filler is perfect to get that beautiful glow. Dermal fillers have quickly become popular over recent years, with K ylie Jenner's perfect jawline and pout adding to the trend. The rise of social media means we want to be photo ready at all opportunities adding to the popularity of the industry. With this continued growth, the safety of procedures is paramount, as we learn more about safe dermal filler. The best dermal fillers are subtle giving a natural look that enhances your beauty and gives you that extra magic. Our aim is transform each look using the highest quality and the safest dermal fillers available in the market.

What makes our dermal fillers stand out?

E.p.t.q uses a high concentration of hyaluronic acid (HA), a sugar found in our skin that holds water and helps it to keep hydrated and plumped. When injected it acts like a sponge to attract water to the skin and plumps the area. HA also boosts the skin's natural collagen meaning fewer touch ups may be needed. Further to this, HA is recognised as safer than other types of dermal filler such as Collagen.

The elasticity of the hive structure in e.p.t.q allows the skin tissue to lift. Our dermal filler uses a very low level of endotoxin minimising the side effects. To reduce pain and discomfort after the procedure e.p.t.q is made to the same PH as the body. Injection pressure is optimised for exquisite procedures to enhance the satisfaction. Our dermal filler uses only raw materials to enhance the safety and quality of the filler. To ensure all products are sterilised e.p.t.q is only produced in an aseptic environment.

How long do e.p.t.q. fillers last?

All fillers are temporary which gives us the positive outlook that the changes we make won't be there forever. However, this means we need to ask how to maintain the look for a longer period of time?

E.P.T.Q fillers last from 12-24 months. Just like any other procedure the results may vary in each independent case. Other factors that can affect how long the filler lasts includes where the filler has been used, how much filler has been used and the speed your body metabolises the filler. We would recommen touch-ups when you start to see your filler decrease in volume.

### E.P.T.Q. dermal fillers

Now you know how our filler works you'll need to know which product you need. We have three types of filler; each achieves

different looks but still all the safest dermal filler for your clients.

S100 - Use it for: correction of the forehead, earlobes, crows feet, temples and ear trophs.

 ${\tt S300 - Use\ it\ for: Lip\ contouring,\ lip\ volume,\ glabella\ lines,\ forehead\ lines,\ cheek,\ chin\ and\ nose.}$ 

S500 - Use it for the correction of age-related changes in skin relief, correction of facial areas with volume deficiency.



E.P.T.Q. S 100

Product ID: A0973

Syringe 1 x 1,1 ml

E.p.t.q. S 100 is a highly purified hyaluronic acid filler with a very low density optimizing the injection pressure to ensure a smoother



E.P.T.Q. S 300

Product ID: A0974

Syringe 1 x 1,1 ml

E.p.t.q. S 300 is perfect for the treatment of middle-range to deep wrinkles in the forehead and correct deep irregularities of the



NO

0**286**19931162089 info@tiar

info@tiantianbike.com

the pain for the patient. The filler is ideal for treating fine lines in the facial area like crow's feet around the eye area, fine wrinkles around the mouth, small wrinkles around the nose and mouth, and tear troughs and temporal hollows. E.p.t.q. S 100 filler restores the skin's natural hydration.

treatment. E.p.t.q. S 100 contains minimizing skin profile, nasolabial folds, perioral wrinkles, and puppet lines. It is ideal to use to shaping and creating more volume to lips. The medium consistency makes it is easy to redefine nose, cheeks, and chin. E.p.t.q.  $\ensuremath{\text{S}}$ 300 is based on highly purified hyaluronic acids.

HUANAN



## E.P.T.Q. S 500

## Product ID: A0975 Syringe 1 x 1,1 ml

E.p.t.q. S 500 has the thickest consistency of the e.p.t.q. line of dermal fillers. This filler is ideal for a more dramatic facial and body sculpturing making it ideal for more highlighted cheekbones and marked jawline. The filler is based on highly purified hyaluronic acid.

SHIJIAZHUANG CITY HEBEI CHINA