



Arasti Plus A Plus White Dermal Filler Hyaluronic Acid Injection 1.1ml

Basic Information

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity:
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms: T/T, Western Union, MoneyGram

Korea

Arasti

White/Plus/A-Plus

CE

1box

box

30USD

5-8 work days

• Supply Ability: 1000box/month



Product Specification

| • Validity: | 3mg/mL 24months |
|--------------------|--|
| Storage Condition: | 2~25°C, And Protect From Direct Sunlight And Freezing |
| • Highlight: | Arasti Plus Hyaluronic Acid Injection, Plus White Hyaluronic Acid Injection, 1.1ml Hyaluronic Acid Injection |



More Images









PRODUCT DESCRIPTION :

ARASTI PLUS

is a premium filler developed by Feel Korea made with purified Hyaluronic acid and infused with lido to provide comfort upon injection. ARASTI PLUS dermal filler is carefully designed to treat medium wrinkles in the nasolabial folds, laugh lines, forehead, nose, and cheek. Once injected, ARASTI PLUS travels into the mid to deep dermis or subcutaneous layer of the skin, effectively enhancing specific parts of the face like the chin and defining the facial features into a more youthful appearance.

ARASTI PLUS is manufactured with good cohesion and elasticity using precise cross-linking technology. The durable structure of this hyaluronic acid dermal filler produces good skin adhesion and natural volume in wrinkled and sagging skin.

It is economically packaged with two pre-filled syringes for your convenience and double the effect at the most affordable price. FORMULATION

Hyaluronic Acid 24mg/mL

Lido 3mg/mL

ARASTI A-PLUS

is the upcoming A-lister in the dermal filler industry with its premium characteristics, making it a star compared to other brands. It is developed by Feel Korea with purified Hyaluronic acid and infused with lido to provide comfort upon injection. ARASTI A-PLUS dermal filler is carefully designed to effectively treat lines and wrinkles found in the forehead, nose, and cheek. It seeps into the deep subcutaneous layer of the skin, which can also be used for chin augmentation.

ARASTI A-PLUS is manufactured through a precise cross-linking technology resulting in a cohesive and stable hyaluronic acid dermal filler. This durable structure can replenish the volume lost in wrinkled and sagging skin.

It is economically packaged with two pre-filled syringes for your convenience and double the effect at the most affordable price. FORMULATION

Hyaluronic Acid 24mg/mL

Lido 3mg/mL

Specification :

| Product name | ARASTI A-PLUS /plus /white Filler | |
|--------------|---|--|
| Application | Fine ,deep wrinkle , Facial Augmentation/Contouring | |
| size | 1.1mL syringe * 2 sterilized needles | |

ARASTI A-Plus



| Component | Cross-linked hyaluronic acid 24mg/mL |
|---------------------------|--|
| Rate of cross- linking | *** |
| Lidocaine | 3mg/mL |
| Validity | 24months |
| Storage condition | $2{\sim}25^{4}{\rm C}_{\rm r}$ and protect from direct sunlight and freezing |
| Packing unit | 1.1mL syringe * 2 sterilized needles(25Gx2) / Box |
| Application part | Deep wrinkles, Facial Augmentation/Contouring |

ARASTI Plus



| Component | Cross-linked hyaluronic acid 24mg/mL |
|-----------------------|--|
| Rate of cross-linking | ** |
| Lidocaine | 3mg/ml. |
| Validity | 24months |
| Storage condition | $2{\sim}25^{\rm s}{\rm C}$, and protect from direct sunlight and freezing |
| Packing unit | 1.1mL syringe * 2 sterilized needles(27Gx2) / Box |
| Application part | Deep wrinkles, Fine wrinkles |

ARASTI White



| Component | Cross-linked hyaluronic acid 24mg/mL |
|-----------------------|--|
| Rate of cross-linking | * |
| Lidocaine | 3mg/mL |
| Validity | 24months |
| Storage condition | $225^\circ\!\mathbb{C},$ and protect from direct sunlight and freezing |
| Packing unit | 1.1mL syringe * 2 sterilized needles(30Gx2) / Bax |
| Application part | Fine wrinkles |

shijiazhuang tiantian technology co., ltd

008619931162089 info@tiantianbike.com info@tiantianbike.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA