shijiazhuang tiantian technology co...ltd

HYARON Dongkook 10syringes prefilled inj Sodium Hyaluronate HYARON Dongkook

Basic Information

Place of Origin: KoreaBrand Name: hyaronCertification: CE

Model Number: 3mlx10pcs
Minimum Order Quantity: 1box
Price: 10USD
Packaging Details: box

• Delivery Time: 1-3 work days

• Payment Terms: D/P, D/A, T/T, Western Union, MoneyGram

• Supply Ability: 1000box/month



Product Specification

Product Name: Hyaron Inj

Material: Hyaluronic Acid

• Content: Sodium Hyaluronate 25mg/2.5ml

Packing Size: 3ml /boxNeedle: 30g

• Function: Remove Wrinkle,anti-aging

• Usage: Skin Care Treatments, Skin Filler/Face Filler



More Images







HYARON prefilled inj Sodium Hyaluronate Dongkook

[Specifications/Box] 10 bottles

[Skin dosage] 3.0ml/bottle, 2 bottles at a time on the face

[Use cycle] once a month, 3 times as a course of treatment

[Maintenance time] About 80 days for the face

[Suitable for skin] People with dull, rough skin, large pores, dry skin, and spots on the face.

[Brand Ring] The third generation of new upgraded version

[Three-generation improvement] HA content> 50 + special whitening ingredients + trace minerals

Introduction to Product Efficacy:

Deep hydration: It continuously replenishes hyaluronicacid in the deep layers of the skin, effectively replenishing the moisture needed by the cheeks, neck, chest, back of the hands and other parts, and enhancing the vitality of the skin.

Restore elasticity: repair the skin barr-ier, promote the regeneration of collagen, gradually restore the elasticity of the skin, and lift and tighten the skin.

Reduce pores: adjust the balance of water and oil, accelerate the speed of skin self-repair, and reduce pores.

Improve skin tone: scavenging free radicals, anti-oxidation, stimulating cell metabolism, quickly expelling melanin from the body, preventing pigmentation, and improving dark yellow and dry skin.

Fade fine lines: Promote collagen proliferation, reduce facial fine lines, stretch crow's feet, Sichuan characters, etc., make the skin smooth and plump.

Application Areas:

Dry skin

First signs of aging

Wrinkles

For alowing effect

Nasolabial folds

Specification

opeomeater :	
Product name	Hyaron inj.
Main Ingredients	Sodium Hyaluronate
packing size	(2.5ml x 10 Syringes)
content	25mg/2.5ml
gauge needle	30g

How to Use:

Hyaron is injected to a depth of 1 mm into the surface layers of the epidermis.

Hyaron gel is administered using a thin gauge needle (30G 1/2").

Usually the procedure is carried out in two stages. The first stage is a course of 4-6 procedures, which are carried out every 10 days. One procedure per month and a half is enough.

Anesthetic cream. The drug administration takes on average 40-60 minutes.

Main Ingredients: Sodium Hyaluronate 25mg/2.5ml

Packing Units: 2.5ml x 10 syringes/box

Hyaron is an excellent option for those seeking an anti-aging solution that delivers visible results quickly. As a mesotherapy product, it can effectively reduce the appearance of wrinkles and nasolabial folds, much like traditional dermal fillers. The product's primary active ingredient, sodium hyaluronate, provides the skin with much-needed hydration and moisture, while also restoring volume and smoothing away wrinkles. Hyaron is a safe and reliable injectable product that can help individuals achieve a more youthful and vibrant appearance

Hyaron is a skin booster for intensive skin .oothing and moisturizing with a formula of highly purified Hyaluronic Acid.Hyaron skin booster is one of the best mesotherapy product which also can be used as a dermal filler. It is a sodium hyaluronate based injectable product with a smooth gel-like consistency that makes it flow excellently. This product is very gentle, making ideal for use both blemish - prone and sensitive skin . The mesotherapy treatment helps rejuvenate and hydrate the skin cells and restore volume on our skin given it a fuller appearance. It so very active in removing wrinkles, boosting skin condition with glowing looks. Being one of the best mesotherapy treatment among others globerly, the effect of the product last long on the skin increasing the time as the treatthmemt goes.

- * It boost the skin appearance by making it tone by reducing pigmentaion .t
- *The product treatment helps eradicate unwanted skin conditions like linear atrophy and acne.
- *The resultant mostly appears in a short period, between 2 to 3 days, and last for a long time
- *Good to go for people looking for dermal fillers with long lasting effect and zero side effect.
- *It improves facial volume, eradicate wrinkles and other spot on the face .















hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA