



Audrey H M 10ml Buttock Enlargement Lip Filler Audrey

Our Product Introduction

for more products please visit us on hadermal-filler.com

Basic Information

- Place of Origin: Korea
- Brand Name: filler
- Certification: CE
- Model Number: 1.1ml 10ml
- Minimum Order Quantity: 1box
- Price: 55
- Packaging Details: box
- Delivery Time: 1-3days
- Payment Terms: T/T, Western Union, MoneyGram, credit card
- Supply Ability: 50000000box/month



Product Specification

- Product Name: Filler 10ml
- Type: 10 MI X1pcs
- Packing: 1 X 1.1 MI Syringe/box
- Typically Performed: 15-30 Mins
- Last Time: 6-18 Months
- Lidoos: 0.3%
- Applications: Facial Lips Nose Check
- Needles: 27G
- Type 1: LYFT
- Delivery Time: 1-2days



More Images



Product Description



shijiazhuang tiantian technology co.,ltd

M Audrey



008619931162089



info@tiantianbike.comhadermal-filler.com

product	dermal filler hyaluronic acid lips facial
Type	1.1 10
size	1.1ml x 1
effection	dermal filler
how to use	injection
usage	lips cheek facial nose filler
usages	dermal filler injection
MOQ	1 piece
eccects	Anti aging wrinkles

AUDREY FILLERS are manufactured by Ostin Bio, Co., Ltd using hyaluronic acid derived from non-animal origins with high purity; hence it is biocompatible when injected into the body. AUDREY FILLERS are composed of highly purified and internationally recognized derivatives of hyaluronic acid that constantly maintain elasticity in the skin.

AUDREY FILLERS are developed with MSFT or Multiple Swelling Formulation Technology, a core method patented by Ostin Bio. It is a method where the hyaluronic acid in the dermal filler uses minimal cross-linking agents. It uses a blend of gels with different swelling ratios to increase hydration and safely maintain volume once injected into the skin. This innovative procedure gives AUDREY FILLERS an excellent volume effect, outstanding tissue adhesion, and high cohesiveness. It also minimizes free radical formation by significantly reducing proteins that may cause hypersensitivity reactions. Upon clinical trials, AUDREY FILLERS shows low unit levels of endotoxin, thus holding an outstanding safety profile.

AUDREY FILLERS come in different variants according to their consistency. The hard variation of the AUDREY Filler is manufactured with a higher viscosity and more volume than the soft type, thus resulting in a more robust shaping capacity after injection. It somehow prevents deformation and can be effective for a longer duration.





Description

1ml with s an efficient dermal filler for correcting moderate to severe lines and wrinkles as well as plumping and shaping your patient's lips. Product is a sterile, transparent gel of stabilized hyaluronic acid of non-animal origin with the addition of 0.3% hydrochloride. 1ml with is supplied in a glass syringe with a Luer lock fitting. Only authorized aesthetic professional in accordance with local regulations should administer with injections.

The packaging content:

- 1 prefilled syringe 1,0 ml
- 2 x needles
- Traceability labels

version is suitable for improvement for fine lines, rinkles and folds. It is highly recommended to use for lips enhancement. The brand is made using patented Non-Animal Stabilized Hyaluronic Acid (NASHA) Technology, which helps retain the natural way that hyaluronic acid molecules connect to each other. This results in a product that produces much longer-lasting results. As well, this filler comes with added , an anesthetic for better patient experience during injection.

Benefits:

Model NO.
radiesse
Efficacy
Wrinkle Removal & Skin Beauty
Delivery Time
3-5day
Transport Package
Box
Specification
1ML
Trademark
revolax
Origin
Korea
Production Capacity
1000000



injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject Revolax with into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject Revolax with into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.injectable gel that contains highly purified hyaluronic acid and to reduce needle pain. Inject into the middle dermis or lip for facial augmentation, wrinkle correction, or lip enhancement. It is manufactured by Across, a leading pharmaceutical company headquartered in South Korea.