



TOXTA 100 Unit Hutox Rentox Innotox 50iu 100iu

Our Product Introduction

Basic Information

- Brand Name: Toxta
- Minimum Order Quantity: 1
- Price: USD 60/box
- Packaging Details: box
- Delivery Time: 1-3 days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 5000 pcs/month



Product Specification


- WhatsApp: +86-19931162089
- Highlight: 100 Unit Hutox Rentox Innotox, Hutox Rentox Innotox 50iu



More Images



Product Description



Toxta is a peripheral muscle relaxant that contains purified botulinum toxin type A (in combination with haemagglutinin) and sodium chloride (excipient). The mechanism of action of the drug is based on the ability of botulinum toxin type A to block the transmission of nerve impulses. Intramuscular injection of the toxin disrupts the synthesis of acetylcholine (the neurotransmitter that provides neuromuscular transmission) and leads to temporary paralysis of the facial muscles. Reducing facial expressions helps to relax the face and smooth out wrinkles and prevent the appearance of new ones.

Main benefits:

Effects on wrinkles on the forehead, around the eyes, between the eyebrows: Smooths wrinkles by inhibiting muscle movement and prevents further wrinkles from forming.

Safe use: Toxta 100U is made from an original botulinum toxin strain from the European National Institute and is cGMP compliant. Low protein content leads to a low risk of immune resistance.

High purity for a faster and more consistent result: High purity (99.8%) and high biological activity ensure fast and consistent results.

Application area:

Mimic wrinkles on the forehead, between the eyebrows and the nasolabial triangle

Eyelid spasm (involuntary contractions of the circular muscle in the eye)

Horizontal and vertical neck wrinkles

Correction of the lip contour (botulinum toxin injections smooth out the curvature of the lips)

bruxism

treatment of hyperhidrosis

Packaging: The drug is produced in the form of a lyophilisate (white powder, which is used to prepare a solution for injection), in glass, tightly closed vials with a volume of 100 units. It should be diluted with 0.9% saline solution.



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA