



## Sculptra Dermal Filler Plla Hyaluronic Acid Dermal Filler

### Our Product Introduction

#### Basic Information

- Brand Name: Aaallright HS550
- Model Number: 1ml
- Minimum Order Quantity: 1
- Packaging Details: box
- Delivery Time: 3-5days
- Payment Terms: L/C, D/P, D/A, T/T, Western Union, MoneyGram
- Supply Ability: 1000box 7days



#### Product Specification

- Highlight: dermal filler plla, sculptra dermal filler plla, sculptra dermal filler

## Product Description

L-poly lactic acid (Sculptra) is a polymer of lactic acid. Synthetic is biodegradable and does not react with our body immunologically, so it is biocompatible. The preparation restores the youthful shape of the cheeks, fills the nasolabial folds, applied to the area of the so-called "Hamsters" or a drooping chin causes flabby areas to be filled and the skin stretched. The preparation is introduced subcutaneously or deeply intradermally in areas where the fat lining is atrophied and lost in volume. The final effect is achieved by rebuilding the collagen structures in the patient's skin, and not by filling. Unlike other treatments that use fillers to temporarily correct imperfections, Sculptra stimulates the skin to produce collagen on its own, which allows for a gradually noticeable and fully natural effect of improving the condition of the skin and rejuvenating the face. Unlike other methods, the obtained effect lasts much longer - up to 2-3 years.



selena wang  
WhatsApp 联系人



### INDICATIONS:

Comprehensive correction of wrinkles  
Rebuilding the volume  
Correction of the face oval  
Loss of skin firmness  
Orbital area defects  
Hand rejuvenation  
Treatment of lipoatrophy - local or general loss of adipose tissue

### APPLICATION:

Face - cheeks, zygomatic bone, nasolabial folds, perioral concavity, marionette folds and chin wrinkles  
Effective for a neck lift  
Can be used on other parts of the body, especially where there is a problem of flaccid skin (e.g. after losing kilograms or after pregnancy)

### TREATMENT PROGRAM:

Usually, full correction requires 2-4 treatments with an interval of about 6 weeks, but it happens that a one-time treatment gives a good result. Durability of results is an individual matter, which depends on the advancement of the aging process, genetic predisposition and lifestyle. The effect lasts an average of two years.

### THE PACKAGE CONTAINS:

Sculptra dry powder is supplied sterilized by sterile filtration in a sterile, oblong clear glass vial with an aluminum cap flush with a rubber stopper and a removable cap. Sculptra powder is prepared with 5 ml of sterile water for injection into a sterile, pyrogen-free suspension.

### EACH VIAL CONTAINS:

150 mg of L-poly lactic acid  
90 mg of sodium carboxymethyl cellulose  
127.5 mg of non-pyrogenic mannitol

Product for professional use only.

By deciding to purchase, you declare that you are a doctor or cosmetologist trained in aesthetic medicine.



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA