



## D plus Elasty Hyaluronic Acid Filler injection

### Our Product Introduction

#### Basic Information

- Place of Origin: Korea
- Brand Name: Elasty
- Certification: CE
- Model Number: D G F
- Minimum Order Quantity: 1box
- Price: 38USD
- Packaging Details: box
- Delivery Time: 5-8 work days
- Payment Terms: T/T, Western Union, MoneyGram paypal ,credit card
- Supply Ability: 1000box/month



#### Product Specification

- Product Name: Elasty
- Type 1: ELASTY Plus D
- Type 2: ELASTY Plus F
- Type 3: ELASTY Plus G
- Content: HA 24mg/mL
- Lido: 0.3%
- Packing Size: 1x1ml
- Whatsapp: 008613315976231
- Highlight: **D plus Elasty hyaluronic acid filler ,  
hyaluronic acid dermal filler injection ,  
D plus Elasty filler with hyaluronic acid**

## Product Description

### Product Description

#### ELASTY

#### HA 24mg/mL

Elasty is a dermal filler contains hyaluronic acid (HA) at a concentration of 24mg/ml, with the addition of lidocaine. Elasty is known for its excellent cohesion and high elastic modulus.

ELASTY Plus with Lidocaine hyaluronic acid dermal fillers are top-of-the-line thanks to their patented Preserved Natural Entanglement Technology (PNET), providing enhanced high elastic modulus and excellent cohesion. These fillers boast high elasticity, high cohesion, and high convenience.

With minimal modification, this dermal filler is highly pure and preferred by consumers. Its high viscosity guarantees maximum satisfaction, making it a perfect choice for quality-focused individuals. High viscosity ensures a proper "fill," achieving optimal results for patients. Additionally, the inclusion of lidocaine minimizes discomfort during the procedure.

**D plus is a dermal filler with medium viscosity. It efficiently fills in wrinkles and folds on the skin and works well on the deep dermal layers. It has been manufactured by cross-linking 3D technology and its high volumetric properties enable resistance to biodegradation.**



ELASTY hyaluronic acid dermal fillers are top of the line in its class with its patented PNET (Preserved Natural Entanglement Technology) which contributes to enhanced high elastic modulus and excellent cohesion. ELASTY touts high elasticity, high cohesion, and high convenience.

With minimal modification, this dermal filler exudes utmost purity when used by any consumer. It makes this filler highly preferred by all consumers. Its high viscosity guarantees maximum satisfaction by any consumer and patient. High viscosity makes for a perfect dermal filler when quality comes into mind. It is because high viscosity guarantees a proper "fill."

ELASTY F has been manufactured by cross-linking 3D technology and its high volumetric properties enable resistance to biodegradation.

### Strengths of ELASTY F Dermal Filler:

- High stage purification process of hyaluronic acid
- Cross-linked fibers ensure a lasting effect
- Painless procedures
- Seamless and no-frills injections
- Easy to bring to the dermatologist

Long biodegradation period (results can last up to a year)  
Ideal for thin skin and sensitive areas  
Low cross-linking technology  
Thin layer washing technology  
High visco-elasticity  
Little post-swelling

## **Product indications:**

ELASTY F filler allows filling in fine wrinkles around eyes, lips, and face wrinkles. It gently enlarges natural volume of tissues.

## **Product specifications:**

### **Hyaluronic acid:**

Contains 24mg/mL

### **Gel type:**

Monophasic gel composition

### **Injection force:**

10~20 N is needed

### **BDDE:**

>0.2ppm

### **Contents:**

1mL x 2 syringes

### **Needle specifications:**

### **Gauge (UTW):**

27G/30G

### **Quantity:**

1EA/1EA

### **Composition (main ingredients)**

HA 24 mg/ml

### **Duration**

The effect of the procedure lasts from 8 to 12 months.

## CASE STUDIES — ELASTY GRAND PLUS



BEFORE



6mL for cheeks, jawline, chin

6mL = 3 boxes Elasty Grand  
Note filler along jawline tightens jowls  
and lifts nasolabial folds



BEFORE



2mL for cheeks & chin

2mL = 1 box Elasty Grand



6mL for cheeks & jawline

6mL = 3 boxes Elasty Grand  
note filler creates neck lift



BEFORE



6mL Elasty Grand + Revolax Fine

6mL = 3 boxes Elasty Grand  
Elasty Grand to enhance cheeks, define jawline, restore chin  
Revolax FINE for under eye circles and bottom lip.  
Note how filler lifts double-chin and sagging around mouth & nose.



BEFORE



1mL for lip augmentation

1/2 syringe (.5mL) used for  
upper lip and .5mL for lower lip



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA