

Ellanse Smart Fillers Next Generation Biocompatible Collagen Stimulator Pcl 30% CMC 70% -C

Our Product Introduction

Basic Information

- Place of Origin: Netherlands
- Brand Name: Ellanse
- Certification: ce
- Model Number: 1ml x 2 per box
- Minimum Order Quantity: 1 piece
- Price: usd 550
- Packaging Details: one box 1 pcs
- Delivery Time: 3-5day
- Payment Terms: T/T, Western Union, MoneyGram,credit card, paypal
- Supply Ability: 1000000pcs/day



Product Specification

- Product Name: Ellanse S/M/L/E Cosmetic Filler Collagen
- Function: Breast Enhancement
- Material: Cross-linked HA, Hyaluronic Acid Filler,
- Application: Breast And Buttock Enlargement, Body Filler,
- Product Name: Ellanse Filler Collagen
- Packing: 1ml X 2 Piece Per Box



More Images



Our Product Introduction

for more products please visit us on hadermal-filler.com

Product Description



shijiazhuang tiantian technology co.,ltd

Ellanse Smart Fillers Next Generation Biocompatible Collagen Stimulator Pcl 30% CMC 70% -C

008619931162089

info@tiantianbike.com

hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA



Ellanse S/M Cosmetic Pcl Filler CMC Treatment Collagen Stimulator Nasolabial Folds Skin Ellanse Collagen Stimulating Filler

Ellansé M 2 x 1 ml Collagen Booster
Ellanse M is a PCL skin filler that can make skin more youthful and radiant, stimulate collagen production and slow down the appearance and formation of shallow lines and wrinkles. It is officially categorized as a Class III medical device and it is usually injected into the sub-dermis. Patient satisfaction rate for this superb cosmetic rejuvenation product is incredibly high and that has a lot to do with its one-of-a-kind properties.

Polycaprolactone (PCL) is a biodegradable substance that makes up almost one-third of the composition of this highly advanced skin rejuvenation product. That is why Ellanse M stands out with its incredible elasticity, viscosity, safety, and durability. Just one treatment with this product generates results that remain active for up to 2 years.

How does Ellanse M stimulate collagen production?

The most unique property of this top-range injection filler is the fact that it increases collagen production within the skin. This reaction is triggered by the PCL microspheres in the dermal implant. The human organism sees them as a foreign body and it replenishes collagen production.

At the same time, however, PCL is biologically compatible with the human body, making the treatment very safe.

The filler consists of a unique patented composition:

70% Carboxymethyl Cellulose (CMC) based gel carrier

30% PolyCaproLactone (PCL) microspheres

Indications:

Improvement of the collagen synthesis

Recovering the shine of the skin

Prevent the appearance of fine lines and

wrinkles

Injection depth:

Indicated for deep dermal and subdermal

implantation

Areas:

Marionette lines

Temples

Cheeks

Nose

Hands

Forehead

Lower face

Brow area

Chin

How does the treatment work?

The treatment can be done in 15 minutes but it can also take longer, depending on the area to be treated.

The professional in charge of your treatment injects the product using an ultra-thin needle under the surface of the skin. For fine lines, the product is injected at the surface of the skin. For marked wrinkles, the injection is made at a deeper level.

Please note that the injection of Juvéderm Ultra requires a medical gesture that can only be carried out by a doctor (doctor specialised in dermoesthetics, dermatologist or plastic surgeon).

