



Handheld Meso Gun Injector 3D Smart Water Injection Pen 3D Smart Injection System

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: Mesotherapy Gun
- Certification: CE
- Model Number: 453
- Minimum Order Quantity: 1
- Price: \$99-110
- Packaging Details: Box
- Delivery Time: 1-3 Days
- Payment Terms: T/T, Western Union
- Supply Ability: 800000000



Product Specification

- Product Name: Portable Mesotherapy Gun
- Feature: Wrinkle Remover, Skin Rejuvenation, Blood Vessels Removal
- Mode: Continuous, Dose
- Material: ABS PC
- Power: Built-in Lithium Battery 7.4V, 850mAh
- Injection Setting: 5 Level (Continuous), 6 Level (Dose)
- Weight: Main Body 220g
- Warranty: 1 Year
- Highlight: **Handheld Meso Gun Injector ,
3D Smart Water Injection Pen ,
RF Beauty Machine injector**



More Images



Product Description

Hyaluronic Acid Ampoule Meso Gun Mesotherapy RF Beauty Machine



Product Description :

Water mesotherapy is to inject hyaluronic acid into skin. Hyaluronic acid contains moisture, 200-300 times as rich as one's body. So that it makes the skin smooth, soft, bright and hydrated after the injection. In order to keep the skin hydrated, ladies always choose varied skin care products. Actually, skin care products only hydrate moisture on the cutin. So that over smearing those products cannot deeply provide moisture for the skin. As we know, acne, wrinkle and inelasticity would appear if skin lacks moisture. Water mesotherapy will solves all the problems by reason of lacking moisture, such as sallow skin, color spot, wrinkle... And it will supply moisture for dermis. Therefore water mesotherapy will works well on , moisturizing and renewing inelasticity.

Function:

- Skin rejuvenation , wrinkle removal and pigmentation removal
- Skin tightening, Enhance skin metabolism, anti-aging
- Activation of collagen cells, enhancing skin elasticity and gloss
- Remove facial fine lines and crow's feet, plain word lines ,etc

Advantage:

- Applicated for endogenous and exogenous injection
- Injection technology which can speed up the treatment process
- Easy to operate, no side effects, but significant effect
- Adjustable injection depth, doses, speed and suction intensity.
- Anti-backflow technology applicated for medicine and blood.
- Stimulate the skin for self-healing, promote skin metabolism

Treatment process:

- Anesthetics
- Fresh-keeping film
- Facial cleaning&disinfection
- Put medicant inside the three-way valve before injection

MEETCARE



Remove eye bag
Only 1minutes

Before → After



Before → After



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA

