



shijiazhuang tiantian technology co.,ltd
hadermal-filler.com

HYALURONIC ACID
DERMAL FILLER
INJECTABLE HYALURONIC ACID

PDO Thread Lift

- The latest and most advanced technology in skin tightening and lifting
- It is made of PDO (polydioxanone) material, which is naturally absorbable in the skin and stimulates collagen synthesis

Beauty 7 Color LED Facial Mask Photon Light Skin Rejuvenation Facial Skin Care LED Light Mask With Neck

Our Product Introduction


Basic Information

- Place of Origin: China
- Brand Name: Mask
- Model Number: Led beauty
- Minimum Order Quantity: 1
- Price: USD 30- \$35/ piece
- Packaging Details: Box
- Delivery Time: 1-3 Days
- Payment Terms: T/T, Western Union, Paypal, Card, Bitcoin
- Supply Ability: 10000000 boxes /day



Product Specification

- Highlight: 7 color LED facial mask, LED light mask with neck, photon light skin rejuvenation mask



More Images



for more products please visit us on hadermal-filler.com

Product Description

Beauty 7 Color LED Facial Mask Photon Light Skin Rejuvenation Facial Skin Care LED Light Mask with Neck



7 Colors LED mask can solve stubborn skin problems



Blue
Sterilization
Inhibit inflammation



Cyanine
Promote metabolism



White
Remove Freckle



Purple
Remove acne scars



Red
Whitening
Skin wrinkles



Green
Row of edema
Reduce skin grease



Yellow
Energy to skin

1. 7 LED lights with multi functions;
2. Easy to operate with full touch screen (3 programs included)
3. Big size for both face and neck;
4. Design patent;
5. Smart design to place the handles on the case;
6. Frequency adjustable on the screen to re-adjust the handle.

Functions:

1. Improve Complexion;
2. Remove Color Spots;
3. Tightening Skin ;
4. Activate Collagen;
5. Anti- Aging Skin .

Red light 630nm	It will stimulate the growth of collagen. Collagen is an essential protein that used to repair damaged tissues and to repair the old tissues. It can remove fine line and shrink pores.
Blue light 470nm	Bacteria in acne contains Porphyrins; because of the high affinity with wave length, the blue light is able to kill Porphyrins. In addition, the blue light has calm action which is very effective for hypersensitivity.
Green light 520nm	The green light can balance the color pigment, and reduce fine lines, improve ageing skin, speed up the healing process of the wound, lightening the skin.
Yellow light 590nm	It has special 590nm wavelength, so it can improve the exchange of oxygen in the cells, to the skin cells to replenish their energy, decomposing pigment, promoting lymphatic drainage, improving the wrinkles, rough skin and etc.
Purple light red+blue	It combines red light and blue light, which has two phototherapy effective . Particularly it has a good effect on repairing the acne,
Clear blue light	It can enhance cell energy gradually, and has a good metabolism facilitating effect.
Laser light (white)	It can penetrate the skin deeply, so to accelerate the tissue metabolism, decomposition color spots, improve the appearance of fine lines and ageing skin.





shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA