

## Hot Cold Hammer Ultrasonic Cryotherapy LED Facial Skin Rejuvenation Lifting Machine

Our Product Introduction

for more products please visit us on [hadermal-filler.com](http://hadermal-filler.com)

### Basic Information

- Place of Origin: China
- Brand Name: Hammer
- Model Number: Hot cold hammer
- Minimum Order Quantity: 1
- Price: USD 35 --\$50/ piece
- Packaging Details: Box
- Delivery Time: 1-3 Days
- Payment Terms: T/T, Western Union, Paypal, Card, Bitcoin
- Supply Ability: 10000000 boxes /day




### Product Specification

- Product: Hot Cold Hammer Ultrasonic Cryotherapy LED Facial Skin Rejuvenation Lifting Machine
- Type: Ultrasonic Cryotherapy
- Effect: Skin Rejuvenation Lifting
- Highlight: **ultrasonic cryotherapy facial machine, LED skin rejuvenation lifting device, hot cold hammer RF beauty tool**



### More Images



#### IMPROVE 5 SKIN PROBLEMS



### Product Description

Hot Cold Hammer Ultrasonic Cryotherapy LED Facial Skin Rejuvenation Lifting Machine

Our F



VISION



## IMPROVE 5 SKIN PROBLEMS



### LARGE PORES

Shrink pores, relieve heat, relieve fatigue, increase skin activity



### SKIN WRINKLES

Ice therapy: regulates cell water content, balances oil secretion, condenses large pores, tightens skin



### EYE BAG EYE DARK CIRCLES

Enhances skin metabolic rate, discharges waste toxins, and has certain effects on eye bags and dark circles

#### 1.Introduction Of Hot Cold Hammer:

##### Cold hammer:

- 1.Be exclusively used for curing a sensitivity skin.
- 2.Can sooner and more availably the opposite department allergic skin have the town calm down function.
- 3.Can immediately remove for five minutes red and swollen, dropsy, stop a itch, go eye bag, solve to bring many skineses unwell face hairdressing sharp weapon because of allergy, let the skin reply a health nature thoroughly.
- 4.For the blush from shame phenomenon in the air condition building, the result is very good.
- 5.Can also used for the assistance treatment that the child gets a fever to have fever.

##### Hot hammer:

- 1.Accelerate product absorption, promote metabolism
- 2.Warm up, make the person feel warm and comfortable, have already relaxed a muscle, equilibrium self-disciplineneve, solve ache,aching
- 3.Stimulate blood circulation acceleration, increase to contain amount of oxygen in the blood, make the pore open, the blood vessel extends to enhance, making the nourishment composition of skin absorb quickly, is the most gentle and comfortable

valid ducting method.

4.Strengthen skin metabolism velocity, make the organization contain hot result to build up, the nature expels discare toxin,having certain result to the cancellation of eye bag, black eye.

**2.Features Of Hot Cold Hammer:**

1. 8 Degree for cold therapy.

2. Portable and easily control Hot surface temperature is up to 42°C Cold surface temperature of the device can be down to 5°C Make face skin tight, smooth and delicated Enhance to help blood circulation, metabolism and relaxation Strengthens collagen;enhances elasticity.

**Functions:**

1.Calming The Skin

2.Shrink Pores

3.Eliminate Swelling

4.Exaggerated Pores


5.Eliminate Bags Under The Eyes


6.Import The Essence

Model NO.	M45
Application	Salon, Home, Beauty SPA
Portable	Portable
Product Name	Hot and Cold Hammer
Function	Skin Rejuvenation Skin Care, Anti-Wrinkles, etc.
Style	Portable
Material	ABS & Stainless Steel
Color	White
Power	50W
Input	AC100-240V 50-60Hz
Use Frequency	15-20min/Time;1day/Time
Hot	Maximum Heating Temperature 40
Treat Area	Face, Eyes, Body
Applications	Beauty Center, Salon, Personal
Shipping	by DHL,UPS, TNT,etc
Delivery Date	2~3 Working Days After Confirming Payment

 **shijiazhuang tiantian technology co.,ltd**

 008619931162089

 info@tiantianbike.com

 haderm:

NO 260 HUANAN ROAD  
SHIJIAZHUANG CITY HEBEI CHINA