



Korean Yvoire Hyaluronic Acid Gel Indicated For Deep Nasolabial Fold Correction

Our Product Introduction

Basic Information

- Place of Origin: Korea
- Brand Name: Yvoire
- Model Number: Contour Classic Volume
- Minimum Order Quantity: 1 piece
- Price: \$20 - \$25
- Packaging Details: box
- Delivery Time: 1-3 Days
- Payment Terms: T/T, Western Union, Paypal, Card, Bitcoin
- Supply Ability: 50000000 pcs/month



Product Specification

- Features: Instant Results
- Ha Concentration: 24mg/ml
- Packing: Prefilled Syringe+color Box
- Suitable For: All Skin Types
- Raw Material: Injection Grade
- Crosslinking: Yes
- Productname: Hyaluronic Acid Wrinkle Fillers
- Duration Of Effect: 6-12 Months
- Skin Type: All Skin Types
- Purpose: Wrinkle Reduction
- Formulation: Gel
- Lid: 3mg/mL
- Private Label: Acceptable
- Shelf Time: 2 Years
- Biocompatibility: High



More Images



Product Description

Korean Yvoire Hyaluronic Acid Gel Indicated For Deep Nasolabial Fold Correction



Product Name	Duration of Effect (Reference)	Core Characteristics & Application Layer	Main Indications & Effects
YVOIRE Classic	Approx. 6–12 months	Small molecule, soft texture; injected into the mid-dermis layer for fine lines and superficial volumization	Improves crow's feet, glabellar lines, and moderate wrinkles; used for lip augmentation and tear trough filling
YVOIRE Volume	Approx. 9–12 months	Medium molecule, good elasticity; injected into the deep dermis for large-area volume restoration and lifting	Fills the cheeks (apple zones) and temples; improves deep nasolabial folds (marionette lines); refines chin contour
YVOIRE Contour	Approx. 12–18 months	Large molecule, strongest support; injected at the periosteum (supraperiosteal) level for structural foundation and skeletal contouring	Facial contour sculpting, such as nose augmentation, chin enhancement, and jawline redefinition; creates sharp, three-dimensional facial features

Product Summary

YVOIRE provides diverse product options meeting needs of doctors and patients.

Product	YVOIRE [®] <i>Hydro</i>	YVOIRE [®] <i>classic</i>		YVOIRE [®] <i>volume</i>		YVOIRE [®] <i>contour</i>	
		<i>classic s</i>	<i>classic plus</i>	<i>volume s</i>	<i>volume plus</i>	<i>contour</i>	<i>contour plus</i>
Purpose	Hydrating	Fine wrinkle correction & Volumizing		Deep wrinkle correction & Volumizing		Volumizing & Contouring	
Product image							
Main particle size	-	400µm		900µm		1400µm	
Composition	Hyaluronic Acid (HA)	Non cross-linked HA		Cross-linked HA		Cross-linked HA	
	Isosane	-	0.3% lidocaine	-	0.3% lidocaine	-	0.3% lidocaine
Enclosed needles	30 G	30 G, 27 G		27 G		21G, 23G, 23G cannula	
Recommended injection layer	Superficial dermis	Mid dermis		Deep dermis		Subcutaneous dermis	
Volume per syringe	1 mL	1 mL	1 mL	1 mL	1 mL	2 mL	1 mL



Why Choose YVOIRE?

YVOIRE offers a complete, tailored range of hyaluronic acid fillers designed to address every layer of facial aging. Its portfolio is built on three key pillars: precision, durability, and safety.

1.Precision by Layer – From the superficial dermis (Classic) to deep structural bone support (Contour), YVOIRE provides a specific product for every indication – fine lines, volume loss, or skeletal contouring – allowing clinicians to treat with artistic control.

2.Long-Lasting Results – With durations ranging from 6 to 18 months depending on the formula, YVOIRE delivers sustained correction, reducing the need for frequent touch-ups and offering excellent value for patients.

3.Proven Safety & Comfort – Manufactured using advanced cross-linking technology, YVOIRE fillers are smooth, cohesive, and designed for easy injection, minimizing swelling and downtime while ensuring natural-looking outcomes.



shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA